

## Atomic Emission Spectrometer AES-3700

### AC/DC ARC Speical Emission Spectrometer

We in domestic take the lead to promote a new combined mode.

Change the photo plate recording model by replacing the photo plate recorder with photo multiplier tube.

Possessing three invention patents, six utility model patents and one analytical software copyright.

2013 BCEIA gold award.

### Features

- ▶▶▶ Combination of AC/DC excitation source, concave grating, PMT detector, make a new combined mode, with an all new optical path, configuration and outlook.
- ▶▶▶ The new light path not only more simple components, but also shorten the optical leagth.
- ▶▶▶ The compact structure arc generation is using a inegraec design, very eary to use.
- ▶▶▶ Can be used to measure geochemical powder samples and non ferrous metallurgy rapidly, accurately and sensivity spcialized analytical method is also afforded.
- ▶▶▶ Automatical setting the exposure time of different element according to the corresponding evaporation curve. The spectral line with different intensity can be measured synchronous.
- ▶▶▶ Specially designed the channel for analytical line,background line internal standard line. Realizing the background. calibration and intemal standard calibration are controlled by program.
- ▶▶▶ The electrode on the excitation device can be imaged on the inspection window which is totally closed, and the electrode position can be adjusted conveniently.
- ▶▶▶ The electrode holders are automatic temperature controlled by water recirculator, only need to adjust up and down, it features a casy of use and long-life.
- ▶▶▶ High integration measure and control system, all the signals can be uniformly controlled. To ensure the real time, efficient and reliable date collection.
- ▶▶▶ The Self-designed analysis software featured with easy-to-use and some special functions, such as subsection integral, background calibration, intemal standard calibration after treatment of analytical data. It also contain powerful database, easy to realizing the data saving, data retrieval and data processing.



## AC/DC ARC Special Emission Spectrometer

### Applications

#### Applications of DC are special emission spectrometer

- ▶▶▶ Measurement of 19 kinds of impurity elements in high purity tungsten oxide.
- ▶▶▶ Measurement of 17 kinds of impurity elements in high purity molybdenum oxide.
- ▶▶▶ Measurement of trace tungsten element in high purity molybdenum oxide.
- ▶▶▶ Measurement of high purity metal powder and metallic oxide, such as copper, aluminium, nickel et al.

#### Applications of AC are special emission spectrometer

- ▶▶▶ International Geochemical Mapping according to IGCP259, IGCP360.
- ▶▶▶ Geochemical Indication for Prospecting.
- ▶▶▶ Stream Sediment and Soil

### Specifications

#### Excitation Light Source

- ▶▶▶ AC/DC are generation
- ▶▶▶ Current 2A~20A
- ▶▶▶ Ignition way: high voltage pulse

#### Excitation Light Source

- ▶▶▶ Monochromator: Paschen-Range Type
- ▶▶▶ Concave grating: radius of curvature 750mm
- ▶▶▶ The reciprocal of dispersion rate: 0.55nm/mm (first order)
- ▶▶▶ Grating: 24001/mm
- ▶▶▶ Wavelength range: 200~500nm
- ▶▶▶ Water-cooling electrode holder (water chiller):  
Temperature: 20°C±2°C
- ▶▶▶ Constant temperature system of monochromator:  
Temperature: 35°C±0.1°C

#### Measure and Control System

- ▶▶▶ The smart measure and control system have multi model designing and multi data collecting function. All the signals can be uniformly controlled.
- ▶▶▶ Measuring style: subsection integral
- ▶▶▶ Reproducibility: RSD≤0.2%
- ▶▶▶ Photomultiplier:  
Voltage: -1000V  
Stability: Better than 0.5% within 8 hours

#### Data controlling System

- ▶▶▶ The Self-designed analysis software featured with easy-to-use and some special functions, such as subsection integral, background calibration internal standard calibration, after treatment of analytical data. It also contains powerful database, easy to realize the data saving, data retrieval and data processing.

#### Dimensions

- ▶▶▶ 1420mm(L)x880mm(W)x1160mm(H)

### Weight

▶▶ 300Kg

### Environment Requirement

▶▶ Environment: temperature:15-25°C  
Relative humidity: ≤75%

## DETERMINED ELEMENTS AND MEASUREMENT RANGE

### Determined elements and measurement range of AC are

- ▶▶ Mainly measured elements
- ▶▶ Determination of trace elements, such as Ag, Sn, Mo, B, Pb, in geochemical samples simultaneously. All the working curve of elements have good linearity respectively, the correlation coefficient of all measured elements are within the scope of 0.9946~0.9999, it also have the function of automatic calculate and display the result and deviation of measurements.
- ▶▶ Measurement of 19 kinds of impurity elements in high purity tungsten oxide  
All the working curve of 19 kinds of impurity elements in high purity tungsten oxide have good linearity respectively, the correlation coefficient of all measured elements are within the scope of 0.9960~0.9998.

### Standard Content and Detection Limit of some element ( $\mu\text{g}\cdot\text{g}^{-1}$ )

Element	Standard Content GBW07701-GBW07709	Detection Limit (DL)
Ag	0.034~10	0.012
Sn	0.28~100	0.24
Mo	0.21~100	0.15
B	2.1~1000	1.65
Pb	2.5~1000	1.30

- ▶▶ Determination of elements, such as Ni, Cr, Co, Be, Cu, Zn, in geochemical sample simultaneously. All the working curve of elements have good linearity respectively, the correlation coefficient of all measured elements are within the scope of 0.9928~0.9992.

### Detection Limit of some other elements ( $\mu\text{g}\cdot\text{g}^{-1}$ )

Element	Detection Limit (DL)	Element	Detection Limit (DL)
Ni	1.5	Re	0.2
Co	0.6	Cu	1.0
Cr	5.0	Zn	10
V	15	Ti	50

### Measurement range of 19 kinds of impurity elements in high purity tungsten oxide

Element	Measurement range%	Element	Measurement range%
Ai	0.000~0.01	Mn	0.0002~0.01
As	0.0005~0.02	Mo	0.001~0.05
Bi	0.00005~0.0024	Ni	0.00015~0.01
Ca	0.00004~0.015	Pb	0.00005~0.0024
Cd	0.00005~0.007	Sb	0.00025~0.01
Co	0.00025~0.015	Si	0.0004~0.02
Cr	0.00025~0.015	Sn	0.00005~0.0024
Cu	0.00003~0.007	Ti	0.00025~0.015
Fe	0.0003~0.02	V	0.00025~0.015
Mg	0.00015~0.01		

- ▶▶ Measurement of 17 kinds of impurity elements in high purity molybdenum oxide  
All the working curve of 17 kinds of impurity elements in high purity molybdenum oxide have linearity respectively, the correlation coefficient of all measured elements are within the scope of 0.9910~0.9999.

## Measurement range of 17 kinds of impurity elements in high purity molybdenum oxide

### Determined gold element

- ▶▶ The chemical spectrometry method, which is frequently used is frequently used to measure the trace gold in geochemical sample, features more sensitivity and time saving.
- ▶▶ Using AES-2700 to measure Au, the detection limit can be 0.3ng.g<sup>-1</sup> the working curve have a good linearity, and the correlation coefficient is equal to 0.9955.


Element	Measurement range%	Element	Measurement range%
Ai	0.0002~0.01	Mn	0.0001~0.008
Bi	0.00005~0.0024	Ni	0.0001~0.008
Ca	0.0004~0.12	Pd	0.00005~0.0024
Cd	0.00005~0.004	Sb	0.0002~0.01
Cr	0.00025~0.012	Si	0.0005~0.012
Co	0.00015~0.012	Sn	0.00005~0.0024
Cu	0.00005~0.007	Ti	0.0003~0.012
Fe	0.0003~0.012	V	0.0003~0.012
Mg	0.0002~0.01	W	0.0053~0.196

### Determined elements and measurement range of DC are

- ▶▶ Measurement of trace tungsten element in high purity molybdenum oxide. The measurement of trace tungsten element in high purity molybdenum oxide is concerned as an acknowledged difficulty in analytical chemistry filed, generally the detection limit only can reach 0.05%. So we innovated “The analytical method of trace tungsten element in molybdenum oxide” if the tungstan element is measured among the content of 0.0053%~0.196%, which featured with a good linearity of requal to 0.9963, the detection limit can reach 0.0027%,which is order of magnitude prompt than usual.

### The working curve of elements of AC/DC are

#### The working curve of Sn in geochemical sample using AC are



#### The working curve of Co high purity tungsten oxide using DC are



## HPLC Servicing, Validation, Trainings and Preventive Maintenance :

- HPLC Servicing** :We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.
- Trainings** :We also take up preventive Maintenance to reduce downtime of HPLC's
- AMC's/CMC** :We offer user training both in-House and at customer sites on HPLC principles, operations, troubleshooting.
- Validations** :We have protocols for carrying out periodic Validations as per GLP/GMP/U SFDA norms
- Instruments** :We offer instruments / Renting Services Modules like pumps, detector etc. on Rent.



## About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO :9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

## About Analytical Technologies



UV/VIS Spectro 2060+ Double Beam



PTIR



Gas Chromatograph 3000



Gas Chromatograph 2979 Plus



Flash Chromatograph



Atomic Absorption Spectrophotometer



Liquid Partical Counter



Optical Emission Spectrophotometer



DSC/TGA



NOVA 2020 plus Automated Bio Chemistry Analyzer



HEMA 2020 Hematology Analyzer



Micro Plate Reader/Washer



Water purification system



Total Organic Carbon



Fully Automated CLIA



NOVA Basic Semi -Auto Chemistry Analyzer



PCR/Gradient PCR/ RTPCR



Blood Gas Analyzer



Random access Analyzer for immunoassay Proteins & clinical chemistry



Semen Analyzer



Water purification carbon

## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



**Analytical  
Foundation**

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [info@analyticalfoundation.org](mailto:info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## Reach us @



**HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)**

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W: [www.ais-india.com](http://www.ais-india.com)  
[www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.