

C-3200

Auto Cell Cytometer



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

This product uses trypan blue staining and two-color fluorescence to stain both live and dead nucleated cells. It can be used in routine cell experiments such as cell killing, cell transfection, apoptosis, cell cycle experiments, bright field cell counting, and GFP positive rate detection.

This instrument is for scientific research testing only and is not intended for use as a medical device. Our company assumes no responsibility or liability for any damage or other losses resulting from failure to operate and use the instrument in accordance.

Features:

- Efficient and convenient to use
- Accurate and quantitative test results
- Combination of fluorescence microscopy imaging and statistical analysis
- Comprehensive provision of cell statistics, average diameter, roundness, agglomeration rate and other information Ability to evaluate cell activity.

Technical Parameters:

1. Single sample loading supports 6 samples, using a disposable cell counting chamber.
2. Optical Parameters: 10x standard objective lens, NA value not less than 0.25.
3. Detection Modes:
 - 1) Brightfield mode, trypan blue detection
 - 2) Fluorescence mode, AOPI detection, GFP positivity detection, RFP positivity detection
 - 3) Fully automatic counting and viability analysis of various cells, eliminating human error.
4. HosSystem:
 - 1) All-in-one touchscreen design, 10.1-inch high-definition 1920*1080 color screen (or higher).
 - 2) Intelligent operating system, standalone app mode, multi-touch, and free image zoom for easy observation of individual cells.
5. Detection Concentration Range: 1×10^4 - 5×10^7 cells/ml.
6. Cell Diameter Range 2-200 μ m, with small cell counting function for precise counting. Cells smaller than 10.
7. Sample volume: 20ul/sample, activity test range 0-100%, coefficient of variation CV<5%.
8. Test Speed: Single well, 3 samples, brightfield <8 seconds, fluorescence <35 seconds, fast testing.

9. Brightfield Illumination: White LED illumination, lifespan >30,000 hours.

10. Fluorescence Illumination:

1) Blue fluorescence center band 480nm

2) Green fluorescence center band 540nm; Supports expandable channel 405nm.

11. Imaging Parameters: Autofocus for optimal image acquisition, camera parameters of at least 6 megapixels, professional-grade CMOS color optical imaging.

12. Device Memory: 500GB solid-state memory, upgradeable.

13. Connection Interface: Ethernet port, two USB 3.0 ports.

14. Software Application Functions: The software has built-in standard cell types and standard cell identification methods to achieve standardized test analysis, including bright field counting, trypan blue counting, AOPI counting, GFP positivity rate detection, RFP positivity rate detection, clumping rate, and nucleation rate detection. The software is free for lifetime upgrades.

15. Data Analysis: Real-time statistical analysis of image information, with high accuracy of results, Database format, results can be saved as PDF

Operation Environment:

1. The instrument must be placed on a dry, sturdy, and non-vibratory surface.
2. Protect the instrument from moisture, dust, smoke, and direct sunlight.
3. The operating temperature must be within a constant range of 18°C to 37°C.
4. Prevent dust and other particles from entering the instrument. Store accompanying reagents and consumables, such as cell counting chambers and pipette tips, in a dust-free environment to avoid affecting image quality.

Operational Precautions:

1. The instrument must be connected to a grounded power source.
2. Operators must be trained and operate this device in a well-managed laboratory. When handling samples, it is recommended to wear protective clothing, especially gloves.
3. Handling of reagents, such as trypan blue, must comply with relevant safety regulations.
4. This device can be used with optional counting chambers. Only uncontaminated, clean, and sterile samples ensure optimal operation.
5. After use, disconnect the USB port.
6. Do not operate this device in a low-temperature environment. If the device is accidentally left in a low-temperature environment, bring it to room temperature and allow it to rest for at least 2 hours before use.
7. Do not place this device on a vibrating laboratory table. Ensure the table is stable and level.

Waste Disposal:

Tissue or cell waste, such as pipette tips, must be disposed of in accordance with relevant biohazardous material handling regulations.

This product uses trypan blue staining and two-color fluorescence to stain both live and dead nucleated cells. It can be used in routine cell experiments such as cell killing, cell transfection, apoptosis, cell cycle experiments, bright field cell counting, and GFP positive rate detection.

Precautions:

To ensure safe and reliable operation of the fully automated cell counter (hereinafter referred to as the instrument), please carefully read and operate the instrument in accordance with this instruction manual. This instrument is intended for scientific research testing only and is not intended for use as a medical device. Our company assumes no responsibility or liability for any damage or other losses resulting from failure to operate or use the instrument in accordance with this instruction manual.

Usage Environment:

1. The device should be placed on a dry, sturdy, non-vibrating surface.
2. The device should be protected from moisture, dust, smoke, and direct sunlight.
3. The operating temperature should be between 4°C and 40°C.
4. Dust and dust particles should be prevented from entering the device. The accompanying reagents and consumables, such as cell counting chambers and pipette tips, should be stored in a dust-free environment to avoid affecting imaging

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



LCMS



Automated Prep-Flash
Chromatography system



Maldi TOF



Optima Gas
Chromatograph



Flash
Chromatograph



DAC
Column



GCMS
3068



UHPLC



HPTLC



Ion Chromatograph



Production
HPLC



HPLC



Column



DLS



Water purification
system

▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility

ANALYTICAL FOUNDATION is a Nonprofit Organization (NGO) for the purpose of:



Analytical Foundation


1. Research & Innovation **Scientist's awards & International Conference visit sponsorship / QC Professional Awards:** Quality life is possible by innovation only and the innovation is possible by research only, hence **ANALYTICAL FOUNDATION** is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research /testing/ publication at info@analyticalfoundation.org

2. Improving quality of life by offering **YOGA Training courses, Work shops / Seminars** etc.

3. **ANALYTICAL FOUNDATION** aims to DETOXIFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Naturecure, Health Care, Ara, Events, Camps etc.

▶▶ A portion of profit of **ANALYTICAL GROUP** Company's goes into **ANALYTICAL FOUNDATION** for the aboce noble cause.

▶▶▶ Reach us @



Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice

HPLC Solutions	MultipleLabs	Infinite Multiplelabs LLP	Analytical Bio-Med	Analytical Foundation (Trust)
----------------	--------------	---------------------------	--------------------	-------------------------------

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E :- info@analyticalgroup.net
info@hplctechnologies.com
info@multiplelabs.com

W:- www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide
USA| UK| Australia| Dubai| Malaysia