

# AUS - 3600

## Automatic Urine Analysis System



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

**Analytical Technologies Limited**

An ISO 9001 Certified Company

[www.analyticalgroup.net](http://www.analyticalgroup.net)

## ▶▶ Test items

### Physical examination

Color, turbidity

### Dry chemistry

Dry chemistry:(11/14 joint tests + 2 ratios):BLD  
BIL,KET,GLU,PRO,NIT,LEU,pH,S-G,UBG,ASC,CRE,Ca,M  
CA,ACR,UCa/Cr

### Formed elements analysis

Red blood cell, leucocyte, cast, crystal, epithelial cell, bacteria and other 31 formed elements.

## ▶▶ AI accurate reconigtion



Red blood cell



Leucocyte



Hyaline cast



Cellular cast



Calcium oxalate  
crystal



Fungus

## ▶▶ Test principle

### Formed element analysis

Based on AI DeepCell engine and digital imaging, simulating the ophthalmoscopic recognition of doctors, AUS - 3600 can intelligently identify, classify , and automatically count the formed elements.

### Dry chemistry test

With using reflection spectrophotometry principle to calculate the reflectivity,then identify the amount of biochemical elements in urine

## ▶▶ **Product Features**

### **Automatic operation during the whole process**

"One press" operation can automatically complete urine sediment sampling, switch objective lens, positioning, focus and scan ; By simulating the action of the doctor ,the test of urine dry chemistry, formed elements, urine color and turbidity is fully automatic without manual intervention during the whole process. Three test modes are optional , dry chemistry mode / formed element mode / dry chemistry + formed element analysis modes can be switched and combined freely.

### **Three core technologies**

#### **i) AI DeepCell engine**

Based on the patented AI DeepCell engine, Reetoo AUS - 3600 Automatic Urine Analysis System can get rid of the shackles of conventional image analysis and significantly increase the test accuracy of urinary formed elements to more than 90%. AI DeepCell engine will learn independently when using, so the recognition rate of complex formed elements will improve consistently.

#### **ii) Formed elements recognition technology of mixed intelligent perception and convolutional neural network**

The depth residual network with adaptive learning mechanism is used to obtain the characteristic parameters of formed elements, and the feature pyramid neural network is used to quantify the features in high and low dimensions, so as to further extract and learn the deep and shallow feature details. Finally, combined with classification convolutional neural network and regression convolutional neural network, the urine formed elements are automatically recognized, classified and counted.

#### **iii) Multi-layer dynamic focus scanning**

In order to solve the problem of some cell blurring due to cell layering in the same field of view, AUS - 3600 intelligently simulates the action of doctors rotating the fine focus spiral under the microscope, with the maximum number of three layers taken under each high-power field. At the same time, patented layered recognition technology is used to identify each layer through using a deep learning algorithm to intelligently identify the optimal image, so the recognition accuracy is greatly improved.



AI DeepCell engine



Formed elements recognition technology of mixed intelligent perception and convolutional neural network



Multi-layer dynamic focus scanning

### High-resolution imaging

The 6.3 million pixel high-quality optical imaging system realizes automatic XY axis scanning ,automatic z-axis focusing, and automatic switching of 10X and 40X HD lenses




### Rotary drum-shaped sampling device

The test strip is put into the test paper bin without considering the direction based on rotary drum-shaped sampling device, which intelligently identifies the test strip direction without manual intervention.



### Disposable microscopic slide

The patented microscopic examination slide is formed by the superposition of carrier slide, sandwich slide and cover slide, which can support the simultaneous tests of 5 samples at a time, avoiding cross infection, ensuring biological safety and higher microscopic examination recognition rate.



### Small and compact model

The world's smallest model provides solutions for microscopic inspection and dry chemistry test, which can save your experimental space and improve the detection efficiency of laboratory unit space.

Size: 670×615×530mm

G.W: about 74KG



NO.	Item	Feature
01	Test mode	Three test modes are optional, dry chemistry mode/formed element mode/dry chemistry + formed element analysis mode (integrated analysis mode)
02	Testing speed	Dry chemistry mode:120T/H, formed element mode:100T/H, integrated analysis mode:60T/H
03	Application scope	Physical tests for urine(color, turbidity), formed element analysis(red blood cell, leucocyte, cast, crystal, bacteria, and others, 37 species in total), dry chemistry analysis(11 /14 joint tests)
04	Sample size	A sample rack has 10 sample tube storage positions, and 5 rows of sample racks can be loaded at most at one time
05	Strip hold capacity	≥100 pcs
06	Microscope chamber capacity	40 slides
07	Objective len	Low power lens: 10X; High power lens:40X
08	Low power lens (LP)	9 visual fields can be collected in 10X objective mode, which is mainly used to test large cells such as cast and epithelial cell.
09	High power lens (HP)	In the 40X objective mode, 9 / 20 / 42 visual fields can be collected, and up to 3 tomograms can be taken in each image, which is mainly used to test small cells, like crystal
10	Carryover	Carry-over rate of formed elements ≤ 0.05%
11	Emergency room	Special emergency room is available, and emergency samples can be tested on arrival
12	Sample processing	Automatically rotate the test tube to identify bar code information
13	Equipment parameter	Size: 670×615×530mm, G.W: about 74KG

## HPLC Servicing, Validation, Trainings and Preventive Maintenance :

**HPLC Servicing** :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

**Trainings** :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

**AMC's/CMC** :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

**Validations** :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

**Instruments** :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



## About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

## Our Products & Technologies




UV/VIS  
Spectro 2080+  
Double Beam



Infra FTIR



Optima Gas  
Chromatograph  
3007



Optima Gas  
Chromatograph  
2979 Plus



Flash  
Chromatograph



Atomic Absorption  
Spectrophotometer



Liquid Particle  
Counter



Optical Emission  
Spectrophotometer



DSC/TGA



Semi Auto Bio  
Chemistry Analyzer



HEMA 2062  
Hematology  
Analyzer



Micro Plate  
Reader/Washer



URINOVA 2800  
Urine Analyzer



Total Organic  
Carbon 3800



Fully Automated  
CLIA



NOVA-2100  
Chemistry Analyzer



PCR/Gradient PCR/  
RTPCR



TOC  
Analyzer



Laser Particle  
Size Analyzer



Ion Chromatograph



Water purification  
system



## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



# Analytical Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [Info@analyticalfoundation.org](mailto:Info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.



**HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)**

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W: [www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)  
[www.ais-india.com](http://www.ais-india.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide