

BSAPA:3480

Bet Surface Area & Porosimetry Analyzers



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

Analytical Technologies Limited is a dedicated manufacturer for gas adsorption analyzers, is a high technology company that combines state-of-the-art hardware and professional theories for nano-materials determination. We originate and serve the Weapons System, relies on advantages of local talented human resource and weapons system technology and devotes into scientific instruments development, production and sale. To research and develop fully automatic and intelligent laboratory determination analyzers, provide international standards, high cost performance and state-of-the-art instruments for R&D institutes, universities and manufacturing enterprises are our unchanged mission.

Our analyzers cover:

- *BET surface area analyzers (flowing gas principle & volumetric adsorption principle)
- *Pore size, pore volume, pore size distribution analyzers
- *High pressure and temperature adsorption analyzers
- *Helium pycnometer density analyzers/ open & closed pore cell analyzers
- *Thermal analysis analyzers (TG/TGA/DTA/DSC etc)

BSAPA-3480 series analyzers work for determining BET surface area mesopore & micropore data by utilizing static volumetric adsorption principle. Be designed from the user perspective and with fully automated operation system, user-friendly interface makes it easy to learn, branded accessories assure stability and prolong its life.

Equipped with 4 analysis ports and 4 pretreatment ports. The pretreatment and analysis ports are separated, can work simultaneously and independently. All testing parameters can be done through software, with full automatic operation and high accuracy data.

They are more suitable for university and commercial labs where have many samples need to be test everyday, with a high-throughput and speedy measurement.

- 10 inch embedded PC with Windows system, can control analyzer operation automatically. Analyzer also has RS232 port for connecting bigger external display.
- 2 testing ports in each side.
- 1 Po port in each side.

- 4L stainless steel LN₂ Dewar, be covered with temperature insulainn material,can support over 72 hours working without refilling liquid nitrogen (LN₂).
- LN₂ tray can realize up/ down smoothly.
- Analyzer main power.
- Temperature probe sensor (in the middle of 3 tubes).
- transparent protective doors be printed with operation notices for testing



Features

Vacuum System

- *unique monolithic manifolds system can decrease connecting points apparently, reduce gas leakage and improve ultimate vacuum.
- *modularity design can configure as customer requests,benefits future functions extension and instrument maintenance.
- *the atlas copco brand bipolar vacuum pump,low noise,stable working,completely provent oil-returning;
- *the ultimate vacuum can reach $4 \times 10^{-2} \text{Pa}$ ($3 \times 10^{-4} \text{Torr}$)



micro-welding manifolds



branded mechanical pump

DEGASSER

Surface area and porosity measurement is closely connected with materials external surface area. Besides, the key of gas sorption is the adsorbates can be efficiently attached onto particle surface or be filled into pores, thus, no more important than samples surface purity & cleanness. The purpose for sample pretreatment is to remove atmospheric contaminants on samples' surface and make room for adsorbates. Prior to any measurement the sample must be degassed to remove water and other contaminants before the surface area can be accurately measured. Samples are degassed in a vacuum at high temperatures. The highest temperature possible that will not damage the sample's structure is usually chosen in order to shorten the degassing time. IUPAC recommends that samples be degassed for at least 16 hours to ensure that unwanted vapors and gases are removed from the surface of the sample. Generally, samples that can withstand higher temperatures without structural changes have smaller degassing times.



Up to 4 samples can be pretreated currently with independent temperature controllers;

- Max temperature is up to 400°C, accuracy $\pm 1^\circ\text{C}$;
- Programmable heating process, step is 1-10°C;
- Flow and vacuum degas modes;
- User-defined analysis gases (N₂ or He are more normally used);
- Can realize fully automatically pretreatment once begin;
- PT100 with resolution 0.01°C, analysis software integrated PID theory.

Technology Features:

1. Adopts stainless steel vacuum system which has perfect sealing performance, high vacuum, stable working and long service duration;
2. Speedy heating process, saving time and improving pretreatment efficiency;
3. Easy installation and uninstallation for sample cells;
4. Unique sample antisplash system can avoid samples be vacuumed into manifolds;
5. Modularity inner structure design, convenient for installation, uninstallation and future upgrade;

6. Each sample pretreatment station has its own adjustable evacuation control.

Environmental:

Ambient temperature: 10 to 50°C;

Maximum relative humidity: 90%.

Electrical: Voltage 100 - 240 VAC; Frequency: 50 or 60 Hz.

Physical:

Height 15.7 inches (40 cm)

Width 11.8 inches (30 cm)

Depth 19.7 inches (50 cm)

Weight 33 lbs (20kg)

Guaranteed Data Accuracy


- * branded silicon thin film capacitive pressure transducer, accuracy can reach 0.1% of real reading, better than 0.1% of F.S. (full-scale) Pressure resolution is 0.016 mmHg.
- * different pressure transducers (1000Torr, 1Torr, 0.1 Torr) be equipped, sectional measurement in pressure range can reduce errors in low vacuum.
- * original stepping coolant level probe system can ensure the coolant level unchanged during the whole analysis process, completely eliminate errors caused by dead volume change.
- * pioneered gas outlet and inlet control system can efficiently prevent sample splash in vacuuming and gas input process, guarantee clean manifolds and steady sample mass, avoid zero and liner drifting caused by transducer's macro-change.



4L stainless steel Dewar



branded pneumatic valves



monolithic manifolds system


Control System

- * branded programmable pneumatic valve system has a strong anti-interference ability, is convenient for installation and uninstallation.
- * separated analysis and pretreatment manifolds can prevent foreign matters to contaminate inner pipelines.


Specs	Models	BSAPA:3480S	BSAPA:3480	BSAPA:3480P	BSAPA:3480TP
		BSAPA:3484S	BSAPA:3484	BSAPA:3484P	BSAPA:3484TP
		BSAPA:3486S			
Sample Pretreatment	Ports	4 pretreatment/degassing ports			
	Temp. Range	RT to 400 °C			
	Vacuum Pressure	vacuuming to 0.5 pa			
	Working Mode	fully controlled by software			
Sample Testing	Principle	volumetric adsorption			
	Ports	4 testing ports work alternately(3480S, 3480, 3480P, 3480TP)			
		4 testing ports work parallel(3484S, 3484, 3484P, 3484TP)			
		6 testing ports work parallel (3486S)			
	Ultimate Pressure	4x10 ⁻² Pa			4x10 ⁻⁶ Pa
	Accuracy	surface area repeatability ≤ ± 1.9%			
		pore size repeatability ≤ 0.02nm			
	Range	surface area only 0.0005 m²/g to no upper limit	surfac area: 0.0005 m²/g to no upper limit pore size: 0.35nm – 500nm		
	Partial Pressure	P/P0 accuracy range 5x10 ⁻⁶ – 0.998			
	Gases	N ₂ , Kr, CO ₂ , Ar, purity 99.999%			
Working Mode	fully controlled by software				
Data Report	single–point BET*	✓	✓	✓	✓
	multi–point BET	✓	✓	✓	✓
	Langmuir surface area	✓	✓	✓	✓
	t–plot surface area	✓	✓	✓	✓
	STSA*	✓	✓	✓	✓
	true density	✓	✓	✓	✓
	BJH* mesopore		✓	✓	✓
	CO2 ice water		✓	✓	✓
	DR* & DA*		✓	✓	✓
	t–plot micropore			✓	✓
	MP micropore			✓	✓
	HK* micropore			✓	✓
	C–Y correction			✓	✓
	SF* micropore			✓	✓
	NLDFT*			✓	✓
Hardwares	Pressure Sensors	0–3 Bar/ 2 pc	0–3 Bar/ 2 pc	0–3 Bar/ 1 pc 0–1 torr/ 1 pc	0–3 Bar/ 2 pc 0–1 torr/ 2 pc 0–0.1 torr/ 2 pc
	Sensors Accuracy	0.15% better than full scale			
	Pipelines	integrated vacuum pipelines			stainless steel micro- welding pipelines
	Valves	solenoid valves + O–rings			pneumatic valves + VCR rings
	Vacuum Pump	mechanical pump			mechanical pump+ turbo pump
	LN ₂ Dewar	4L stainless steel Dewar, can work over 72 hours without intervention.			
	Sample Tubes	quartz material			
	Dimension	net: 75x60x90 cm,80 kgs gross: pkg 1: 92x73x110cm,115kgs; pkg 2: 60x75x80cm,70 kgs.			
Extra Accessories*	helium cylinder	purity over 99.99%, volume can be 10L, 20L or 40L cylinder			
	nitrogen cylinder	purity over 99.99%, volume can be 10L, 20L or 40L cylinder			
	liquid nitrogen (LN ₂)	should prepare at least 10L for calibration, it is consumblething.			
	balance (option)	precision should be 0.0001 g			

Remarks:


1. BET – Brunauer, Emmett, and Teller
2. STSA – Statistical Thickness Surface Area
3. BJH – Barrett, Joyner, and Halenda
4. DR – Dubinin–Radushkevich
5. DA – Dubinin–Astakhov
6. HK – Horvath–Kawazoe
7. SF – Saito and Foley
8. NLDFT – Non-local Density Functional Theory
9. extra accessories should be prepared by buyers because they (except balance) cannot be shipped by air.




Single point BET & multi-point BET report



Adsorption & desorption isotherms




BJH (mesopore) adsorption



SF (micropore) pore distribution

Summary Report			
Surface Area			
SN	Items	Description	Result
1	Single point surface area	P/P0=0.2912	284.507898(m ² /g)
2	BET Surface Area	Point range:0.0087 - 0.1085	233.630226(m ² /g)
3	Langmuir Surface Area	Monolayer adsorption model calculation	274.943275(m ² /g)
4	t-Plot Micropore Area	A.Harkins-Jura Adsorbed thickness (nm):0.3638 - 4.4 nm	292.742876(m ² /g)
5	t-Plot External Surface Area	SBET - Sinteric	40.260553(m ² /g)
6	BJH Adsorption cumulative surface area	Pore width (nm): 1.8986 - 305.6535	23.326455(m ² /g)
7	BJH Desorption cumulative surface area	Pore width (nm): 2.0047 - 305.6535	35.222558(m ² /g)
Pore Volume			
SN	Items	Description	Result
1	Single point adsorption total pore volume	P/P0 = 0.9817 - total pore volume of the critical pore width here: 0.0035-0.0035	0.172248(cm ³ /g)
2	t-Plot micropore volume		0.117341(cm ³ /g)
3	SF (S, H, H) (t-Plot micropore volume)	P/P0 = 0.0064 - total pore volume of the critical pore width here: 0.0035-0.0035	0.132159(cm ³ /g)
4	BJH Adsorption cumulative volume	Pore width range (nm): 1.898773 - 305.653454	0.040648(cm ³ /g)
5	BJH Desorption cumulative volume	Pore width range (nm): 2.004720 - 305.653454	0.050699(cm ³ /g)
Pore Size			
SN	Items	Description	Result
1	Total adsorption average pore width	By t-Plot A Adsorption BET specific surface area value	2.066524(nm)
2	BJH Adsorption average pore width	By t-Plot A Adsorption cumulative pore surface area	7.381498(nm)
3	BJH Desorption average pore width	By t-Plot A Desorption cumulative pore surface area	5.314623(nm)
4	BJH Median Pore width	Pore width range (nm): 2.037987 - 305.653454	2.037987(nm)
5	SF Median Pore width	Pore width range (nm): 0.4367 - 1.9902	0.557674(nm)

Comprehensive report summary



t-plot report

APPLICATIONS

Battery Materials:

Cathode materials:

NMC (NCM) – Lithium Nickel Cobalt Manganese Oxide (LiNiCoMnO_2)

LFP – Lithium Iron Phosphate (LiFePO_4/C)

LNMO – Lithium Nickel Manganese Spinel ($\text{LiNi}_{0.5}\text{Mn}_{1.5}\text{O}_4$)

NCA – Lithium Nickel Cobalt Aluminum Oxide (LiNiCoAlO_2)

LMO – Lithium Manganese Oxide (LiMn_2O_4)

LCO – Lithium Cobalt Oxide (LiCoO_2) lithium titanate etc.

Anode materials: typically graphite, others like silicon, tin and its alloys. Separators: typically microporous polypropylene or polyethylene film

Petrochemical Industry:

carbon black, surfactant, white carbon black, pesticide, titanium dioxide, additive, alumina, molecular sieve catalyst, resin, carbon fiber, etc

Medical Implants:

montmorillonite powder, magnesium stearate, micropowder silica gel, talc and other lubricants, diluents such as starch and sugar, aluminum hydroxide and other pharmaceutical excipients

Environmental Protection Materials:

pine, activated carbon, bamboo charcoal, chemical carbon, zeolite molecular sieve, etc

Nano Materials:

MOF materials, nanometallic materials, nano silica, nano zinc oxide, carbon nano tubes, nano aluminum oxide, nano silicon nitride, graphene, etc

Other Powder and Granular Materials:

magnesium oxide, zirconia, calcium oxide, metal powder, ceramics, mineral powder, etc



HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Partical
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide