

GT-3001

Ghost Peak Trap Column



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

Ghost peak is a collective term for chromatographic peaks that appear sporadically or inexplicably during chromatographic analysis. There are many reasons for ghost peaks, such as impurities in the mobile phase, residual samples in the pipeline, other impurities in the mixer, etc., but ghost peaks caused by impurities in the mobile phase are the most common.

Ghost Peak Trap Column Function:


Can effectively inhibit the appearance of ghost peaks in the gradient analysis to ensure the accuracy of experimental results.

Reduces baseline noise and increases detection sensitivity.

Ghost Trap Column can effectively trap impurities in the mobile phase, flow path and mixer, especially can significantly reduce the impact of ghost peaks on analytes in the analytical chromatogram during gradient elution, thus shortening the time for method validation and impurity analysis, improving work efficiency and ensuring the accuracy of results.

Connection of Ghost Trap Column:

The Ghost Trap Column is installed between the pump system and the injection system in the following position.



Use of Ghost Trap Columns:

The ghost trap column should be connected after the convergence of the gradient mixer or pump and before the injector, otherwise it may cause strong adsorption of the sample. Before using the ghost column, please do not connect to the back of the system, rinse with 85% methanol solution, 1mL/min for 30min, and then connect to the equipment. Before and after the use of buffer salt, please pay attention to the use of a high proportion of water transition, to avoid the buffer salt crystals precipitation, resulting in packing blockage. When the trapping effect is found to be less than satisfactory, try to rinse the column with pure acetonitrile.

Information:

Product Name	Specification	Volume	Pressure	Application equipment
GT-3001	4.6 x 50mm	800μL	40MPa	HPLC
Ghost Trap Column	2.1 x 50mm	200μL	100MPa	UHPLC

Practical cases:

Case 1: Trifluoroacetic acid solution

Mobile phase : A: 0.05% TFA aqueous solution


B: 0.045% TFA acetonitrile solution

Column temperature: 40°C .


Wavelength: 275nm

Gradient:

Time	0	10	12	20	22	25
Mobile phase B	30	30	56	65	95	30
Flow rate (mL/min)	1.5	1.5	1.5	1.5	2.0	1.5



Blank spectrum without ghost peak trapping column



Blank spectrum of column with ghost peak

Case 2: Phosphoric acid solution

Mobile phase : A: 0.1% aqueous phosphoric acid solution

B: methanol

Column temperature: 35°C


Flow rate: 1.0mL/min

Wavelength: 225nm


Wavelength: 225nm

Gradient.

Time	0	10	25	30	31	38
Mobile phase B	60	70	90	90	60	60



Blank spectrum without ghost peak trapping column



Blank spectrum of column with ghost peak

Precautions for use of Ghost Trap Columns:

The ghost trap column should be connected to the gradient mixer or pump after the sink outlet, the lagging volume is the same as the column volume.

The ghost trap column must be connected before the injector, otherwise the sample may be adsorbed; if ion-pair reagents are used in the mobile phase, the trap column may have certain adsorption effect on the ion-pair reagents, which may affect the retention time of the components or the peak shape.

Thoroughly flush the system lines with the mobile phase before connecting the column (close to the final concentration in the gradient analysis).

The actual lifetime of the trap column will vary depending on the mobile phase used. Not all impurities in the flow path can be removed, so please use it according to the actual situation.

Case 3: Ammonium phosphate solution

Mobile phase : A: 5mmol/L aqueous diammonium hydrogen phosphate:10mmol/L dihydrogen phosphate Ammonium phosphate aqueous solution=1:1
B: Acetonitrile

Column temperature: 30° C


Flow rate: 1.2mL/min

Flow rate: 1.2mL/min


Wavelength: 275nm

Gradient.

Time	0	30	40	40.1	47
Mobile phase B	25	75	75	25	25



Blank spectrum without ghost peak trapping column



Blank spectrum of column with ghost peak

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



LCMS



Automated Prep-Flash
Chromatography system



Maldi TOF



Optima Gas
Chromatograph



Flash
Chromatograph



DAC
Column



GCMS
3068




UHPLC



HPTLC



Ion Chromatograph



Production
HPLC



HPLC



Column



DLS



Water purification
system

▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility

ANALYTICAL FOUNDATION is a Nonprofit Organization (NGO) for the purpose of:



Analytical Foundation

1. Research & Innovation **Scientist's awards & International Conference visit sponsorship / QC Professional Awards:** Quality life is possible by innovation only and the innovation is possible by research only, hence **ANALYTICAL FOUNDATION** is committed to identify such personallities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research /testing/ publication at info@analyticalfoundation.org

2. Improving quality of life by offering **YOGA Training courses, Work shops / Seminars** etc.

3. **ANALYTICAL FOUNDATION** aims to DETOXIFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Naturecure, Health Care, Ara, Events, Camps etc.

▶▶ A portion of profit of **ANALYTICAL GROUP** Company's goes into **ANALYTICAL FOUNDATION** for the aboce noble cause.

▶▶▶ Reach us @



 **Analytical**®
Technologies Limited

Subsidiaries



Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice

HPLC Solutions

MultipleLabs

Infinite Multiplelabs LLP

Analytical Bio-Med

Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E :- info@analyticalgroup.net
info@hplctechnologies.com
info@multiplelabs.com

W:- www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide
USA| UK| Australia| Dubai| Malaysia