


Infrared Forehead Thermometer



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Instructions

- The temperature measurement range is 95~45 °C. when the temperature is low, "LO" will be displayed. Above 45 °C, "HI" will be displayed
- Open the battery cove, install the 9V battery, and close the battery cover. Press and hold the key to power on, enter the measurement mode after displaying "LO" after 3 seconds. In the measurement mode, each time you press the key, the temperature can be measured once in about 1 second. "SCAN" is displayed on the LCD. After the switch is released, the current measured temperature value is locked on the LCD.
- ALARM (high temperature alarm) function settings: Press the ALARM key, the letter "R" is displayed on the left side of the LCD, use the up and down keys to select the alarm temperature point, and then press the ALARM key to exit. When the measured temperature exceeds the set alarm temperature, the buzzer will issue a "Bi, BI" alarm sound
- The default is no backlight. Press the LIGHT button to turn on the backlight function, and then press the LIGHT button to turn off the backlight function.
- When the button is not operated, it will automatically shut down in about 60 seconds.

►► Specifications

Measuring Temperature Range	35~45 °C
Measurement Accuracy	± 0.3 °C Operating temperature: 23 °C ± 3 °C
Emissivity	0.95 Preset
Response Time	1~2s
Response Wavelength	8-14um
Resolution	0.1 °C
Measuring Distance	3cm ± 0.5cm
Recommended Use Environment	15 - 35 °C
Relative Humidity	10-95 % RH Non-condensing
Storage Temperature	-20~60 °C, Except battery
Weight/Size	~ 130g ; 146 x 80 x 38 mm
Using Electric	9V 6F22 battery
Usage Time	~ 12 Hours

►► Features

- Non-contact high-precision measurement
- can set alarm temperature
- Automatic data retention and automatic shutdown
- Resolution is 0.1 °C
- LCD display with backlight


1 Measurement switch: When the switch is held, the temperature value is displayed, and (SCAN) is displayed at the same time. When the switch is released, it is changed to "HOLD" and the temperature display, and the data is automatically maintained. It will automatically shut down after 20 seconds without operation.

2 Battery cover

3 ALARM : temperature measurement alarm setting key

4 Backlight switch and up adjustment key

5 Display




►► Instruction For Use

1 Open the battery cover, properly install the battery, and turn on the battery switch button to turn on. The screen displays "--", wait 3 seconds to start measurement.

2 Display Symbol:

- a. Data lock symbol
- b. Measuring symbols
- c. Backlight on symbol
- d. Battery level indicator
- e. Temperature Measurement unit symbol
- f. Temperature Measurement unit reading



3 Locating the temperature measurement point: press and hold the switch button, and simultaneously point the measurement at the measurement target to measure.



►► **Details Show**

Three panel hood design

Protective hood design; elastic cap mouth;
Effectively block dust and biological contaminants



Zipper design

Simple and easy to wear and take off



Elastic cuff design

Comfortable and easy to work



Elastic waist design

Elastic waist design, more comfortable to wear,
but also to meet the needs of different staff figure.



Coveralls design

Coveralls design provide additional protection,
breathable and waterproof, no cross infection.



HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Particle
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility



Analytical
Foundation

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org
2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.
3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide