

MP - 3000Series

MERCURY POROSIMETER



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Characterization:

- The knowledge of porosity, pore sizes and pore volume is of fundamental importance for the characterization of porous materials. Mercury porosimetry is the most widely used method for determining the pore size distribution of accessible macro- and mesopores in solids.

The technique is based on the pressure-dependent intrusion of mercury as a non-wetting liquid into a porous material. Using the Washburn equation, the corresponding pore size is calculated from the applied pressure.

As a provider of state-of-the-art solutions in the field of particle characterization and gas adsorption, has now expanded its extensive portfolio with a series of devices specifically for the field of mercury porosimetry. The consists of state-of-the-art instruments for measuring pore size distribution, pore volume, specific pore surface area, density and particle distribution of finely divided and porous materials.

The mercury porosimeters precisely and quickly detect all accessible pores in the size range from 1 millimeter to 3.6 nanometers with stepwise equilibrium-controlled pressure increase from a vacuum up to 414 MPa.




►► Features :

- * Fully automatic vertical filling under constantly high vacuum
- * High resolution enables detection of up to 20,000 measuring points
- * Secure and full functionality without gas connection and liquid nitrogen
- * Low space requirement due to compact design
- * Vertical arrangement of the dilatometers ensures safe handling
- * Efficient re-use of mercury via cleaning set

MP-3001/3002

MP-3003

- 1 Status display
- 2 Manual control
- 3 High-pressure oil supply
- 4 High-pressure autoclave
- 5 Mercury reservoir
- 6 Analysis and vacuum port



►► **Accurate Determination of The Pore Size Distribution:**

MP-3001




MP-3002



MP-3003



Pressure area



The 3000 Series mercury porosimeters from Microtrac MRB reliably and reproducibly measure pore diameters from 1 millimeter to 3.6 nanometers at 414 mega-Pascals. “Pascal” is not only a pressure unit, but also stands for the equilibrium-controlled and optimized handling of the pressure build-up through the so-called “Pressurization by Automatic Speedup and Continuous Adjustment Logic”, or P.A.S.C.A.L. for short. This automatic control is regulated by the real pore system and permits shorter measuring times under guaranteed equilibrium conditions as well as the detection of all pores within the specification - and this with up to 20,000 measuring points per analysis. Since only three types of dilatometers are sufficient for all measurement tasks and neither gases nor liquid nitrogen are required, running costs can be kept significantly low. In addition, the 3000 Series LP low-pressure porosimeter is easy to use and has an extended measuring range up to pore sizes of 1 mm. The vertical degassing and filling with mercury on the 3000 Series LP allows the degassing pressure to be adjusted, making it possible to measure moist samples without changing the moisture content of the material. In this way, even moist concrete samples and solvent-containing porous green bodies can be measured unaltered. The mercury porosimeters from Microtrac MRB are available for different pore ranges:

3000 Series MP-3001 (1,000 - 3.25 μm)

3000 Series MP-3002 (15 - 0.0065 μm)

3000 Series MP-3002 (15 - 0.0036 μm)

►► DILATOMETERS & ACCESSORIES

Available accessories:

- **Dilatometer (sample vessels)**

Dilatometers in different sizes for the MP-3000 series

- **Ultra-Macropore Set (UMP)**

Extends the measuring range of the 3000 Series LP for pore & particle size determination

- **Porosimeter calibration set**


For easy calibration of capacitive volume measurement

- **Mercury cleaning set**

Enables the efficient reuse of mercury



The Professional dilatometer offers an easy & safe handling via screw cap and integrated opening aid




DILATOMETER	Professional	Standard	UMP (Ultra Macropore)
Capillary diameter	3 mm	3 mm 6 mm	6 mm
Sample type	Powder // solids	Powder // solids	Powder, solids
Max. sample size solid state body (d x h)	11 x 25 mm	12 x 46 mm 25 x 25 mm	25 x 25 mm
Compatible with	LP, MP, HP	LP, MP, HP LP, MP	LP
Dilatometer volume	8 cm ³	15 cm ³ 35 cm ³	50 cm ³

►► POWERFUL & INTUITIVE EVALUATION SOFTWARE

Features

- Pore size distribution (differential, integral and as histogram)
- Pore volume, porosity
- Particle size distribution
- Bulk density and apparent density
- Specific surface
- Fractal dimension
- Tortuosity
- Permeability
- Frost resistance of concrete
- Compliant to DIN ISO 15901



The PoreInspect software for the 3000 Series instrument series of Microtrac MRB offers a variety of functions. The software monitors and controls up to 4 measuring instruments connected via LAN completely independently with a PC and allows for individual control of all instrument parameters in real-time.

The data evaluation of the PoreInspect software not only takes the exact blank value correction into account, but also allows for a check and correction of possible sample compressibilities, guaranteeing most reliable results. Extensive evaluation options enable the selection of different pore models and their presentation as graphics and histograms. The calculations are freely selectable depending on the application and include, for example, the representation of fractal dimensions as well as calculations of tortuosity and permeability. Overlay, statistical evaluations, the creation of a method catalogue and data export are further features. The stored raw data are available at any time.

Model	MP-3001	MP-3002	MP-3003
Function	Degassing, Hg filling, low-pressure porosimetry	High-pressure porosimetry	High-pressure porosimetry
Pressure range	Vacuum up to 450 KPa	0.1 - 228 MPa	0.1 - 414 MPa
Resolution	0.001 kPa in measuring range: vacuum to 0.1 kPa 0.01 kPa in measuring range 0.1 - 450 kPa	0.001 MPa in measuring range 0.1 - 100 MPa 0.01 MPa in measuring range 100 - 228 MPa	0.001 MPa in measuring range 0.1 - 100 MPa 0.01 MPa in measuring range 100 - 414 MPa
Pressure detection accuracy	better than 0.1% F.S.	better than 0.1% F.S.	better than 0.1% F.S.
Pore size diameter	180 - 3.25 µm (UMP: 1,000 - 3.8 µm)	15 - 0.0065 µm	15 - 0.0036 µm
Particle size diameter	330 - 15 µm (UMP: 3,000 - 15 µm)	40 - 0.015 µm	40 - 0.01 µm
Max. detectable volume	0.5 cm ³ - 2 cm ³	0.5 cm ³ - 2 cm ³	0.5 cm ³
Volume detection accuracy	better than 1% F.S.	better than 1% F.S.	better than 1% F.S.
Max. number of measuring points	10.000 intrusion 10.000 extrusion	10.000 intrusion 10.000 extrusion	10.000 intrusion 10.000 extrusion
Weight & dimensions (W x D x H)	55 kg (121 lbs) 40 x 67 x 80 cm	68 kg (150 lbs) 40 x 67 x 80 cm	68 kg (150 lbs) 40 x 67 x 80 cm

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph




Atomic Absorption
Spectrophotometer



Liquid Partical
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**[®]
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.