

P-OES-3600

Portable Optical Emission Spectroscopy



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

Summary:

Portable optical emission spectroscopy that can be used in many applications like the metal producing, processing, and recycling industries. P-OES-3600 mobile metal analyzer has the advantages of high analysis accuracy, good analysis stability, convenient operation, wide application, small size, light weight, etc. The use of widely used and trusted optical emission spectrometry(OES) technology to make P-OES3600 mobile metal analyzer results more reliable.

P-OES-3600 is your choice to analyze the contents of nonmetallic elements such as C, P, S, B, Sn, As and other metals when it is inconvenient to cut large metal parts. It can analyze various forms of metal materials, such as pipes, bars, valves, welds, tanks, castings, etc. P-OES 3600 spectrometer can meet your analysis requirements, whether in or out of the factory, on the high ladder in the waste yard, underwater diving tank or chemical plant.

Scope of supply:

No.	Item	Specifications	QTY	Unit	Remarks
1	Mobile Optical Emission Spectrometer	P-OES	1	Set	Include in quotation.
2	Mobile cart		1	Set	
3	Laptop computer		1	Set	
4	Tungsten electrode		1	Pc	
5	Electrode brush	Φ6	2	Pcs	
6	Optical lens	Φ10×2.65	2	Pcs	
7	Degreasing cotton		15	Gram	
8	Lithium Battery	24V/10A	1	Set	
9	M4 Internal hexagonal wrench		1	Set	
10	Spark stand screw	M4	2	Pcs	
11	Argon gas pressure reducing valve		1	Pc	
12	Argon gas pipeline		2	M	
13	Spectrum calibration samples		1	Group	
14	Hand polishing machine		1	Set	
15	Manual		1	Set	

Optional item	Qty	Note
Argon purifier	1 set	When argon gas is not as pure as 99.999%.
Type samples	pieces	Essential for foundry

Analytical Program

		Fe-001	Fe-003
No	Elements	Low Alloy Steel	Cr/Ni Stainless Steel
1	C	0.006-1.3	0.008-2.5
2	Si	0.01-2.9	0.09-4.0
3	Mn	0.03-14	0.12-16
4	P	0.002-0.12	0.003-0.3
5	S	0.002-0.46	0.001-0.4
6	Cr	0.01-12.5	7.4-32
7	Ni	0.004-4.4	0.8-40
8	Mo	0.004-1.76	0.08-4.2
9	Al	0.003-0.5	0.005-1.7
10	Cu	0.002-0.7	0.05-4.5
11	Co	0.001-0.5	0.008-0.62
12	Ti	0.002-0.5	0.005-1.1
13	Nb	0.002-0.53	0.02-2.0
14	V	0.003-0.9	0.02-0.58
15	B	0.006-0.02	0.007-0.02
16	Sn	0.001-0.09	0.003-0.05
17	As	0.001-0.1	0.004-0.04
18	Sb	0.002-0.02	
19	Fe	Reference	Reference

Technical Parameter:

Optical System:

Detector system: high-performance CMOS detector display

Wavelength range 165-585 nm

Spark Unit:

Argon Sparge Electrode The standard configuration length of optical fiber is 3.5m (it can be lengthened according to user's requirements).

Equipped with HPC fiber

Probe is exquisitely designed and easy to clean

UV-Touch Probe:

Used to detect typical ultraviolet elements (C, P, S, B, Sn, As, etc.)
Configure the network information connected to P-OES-3600
Cloud platform operation

Spark Source:

High energy spark light source
Built-in energy electronic regulation and steady-state device
Unidirectional discharge
Computer control of light source parameters
Frequency: 100-1000 Hz
250-300V
Maximum pulse current 80A,

Computer and Readout System

Built-in industrial control host
Battery powered
LAN, USB port open
Built-in HD LCD display
Rugged and durable screen, no keyboard required, built-in mouse, easy to operate

Spectral Analysis Software

Windows 10 operating system
Dedicated spectrometer software for multitasking operation:
Raw Material Sorting
Automatic calibration control is carried out according to the given deviation and excitation times.
Electrode cleaning according to set data
Display and corresponding deviation
Detecting element results storage and printing
Freely editable reports
Display light intensity value, intensity ratio, recalibration and corrected intensity ratio
Analysis result percentage content display
Additive or Multiplicative Interference Correction
Matrix correction
Automatically select spectral lines
Sample identification input
Analysis data can be transferred to an external computer
Connecting external printers and monitors
Analysis data is stored in standard formats such as Excel, ASCII, etc.

Technical Data

Cart dimensions:

Length: 780mm, Width: 610mm, Height: 1050mm

Weight: about 20kg (including excitation tip and box. Argon cylinder and mobile cart are not included)

Power requirements:

Configure 220V power charger Argon gas requirement:

Argon purity 99.999%

Inlet pressure 0.3Mpa

Technical Features

1. Multi-matrix, All-Element Analysis Technology

Designed for rapid identification and accurate analysis of metal materials in industrial field, it can accurately determine C, P, S and As in steel. The matrix of Fe, Cu, Al, Ni, Ti, Sn and Co can be analyzed. It can be comparable to the high-precision analysis of laboratory spectrometers and can also be used for quick material sorting with analysis time of only a few seconds. It can analyze metal materials and products such as pipes, bars, valves, welds, oil tanks and large castings. The compact design makes the equipment multifunctional and can meet different application requirements.

2. Unique Optical Design Technology

The optical system of the equipment is integrated with CMOS detector technology, and all the spectral lines can be read. The compact and lightweight design can truly realize the "portable" and firm mechanical structure, and can adapt to various vibration environments. Very wide spectral line range and very high resolution can be reflected from the number and precision of elements in the test results.

3. Intelligent Calibration Function

Eliminate the spectral drift caused by the change of temperature and pressure; Full-spectrum technology automatically and quickly addresses, which avoids laborious manual operation and human error.

4. Full-Spectrum Direct Reading Technology

The equipment can easily increase spectral lines and matrix. Multi-matrix analysis can be realized without adding hardware facilities.

Customers can increase the types of analytical matrix and elements on the equipment according to the actual needs in the use process, which can be conveniently completed at the user's site.

5.High-Performance Special Optical Fiber Technology

Since the spectrometer is in a high-energy and high-temperature state during the excitation process, the performance of ordinary silica fiber will be drastically reduced due to overexposure under such working conditions. The breakthrough high-performance special fiber technology solves this problem, ensuring that the fiber can work in the best condition for a long time and stably.

Since the optical fiber can maintain long-term stable work during the signal transmission process, the equipment has a good effect in testing the application of ultra-low carbon and other trace elements, and it is convenient to complete the sorting and identification of steel on site without resorting to laboratory equipment.

6.Unique Jet Electrode Technology

In the excited state, a jet of argon gas is formed around the electrodes. The argon gas flow around the excitation point ensures that the excitation process is not disturbed by the outside world, saves the use of argon gas, and reduces the use cost of customers.

7.Multifunctional Adapter

Wire (with a universal adapter, the smallest wire can be analyzed 1 mm) Pipe, bar, valve, weld, etc.

8.Intelligent Spark Design

The design of the spark tip is ingenious, fully considering the needs of practical applications. It can detect low content of C, P, S, B, As and Sn in medium and low alloy steel and stainless and the main functions of the equipment can also be directly controlled.

9.Direct Power Supply from Lithium Battery

Fully consider the actual application situation in the industrial field, and ensure that customers can maximize the working efficiency of the equipment. It is powered by a new lithium-ion battery, with sufficient power, convenient use, safety and reliability. The new lithium-ion battery can stand by for up to 10 hours and can analyze hundreds of samples.

10.Operating Software Based on Win10 System

The operating software is easy to learn, with both Chinese and English versions. Operators without any knowledge of spectroscopic equipment can use the equipment with only simple training.

11.It Embodies the Combination of Humanization and Modern Scientific and Technological Concepts

The design of the equipment includes the latest CMOS and readout system technology, as well as the Win 10 operating system, so that the subsequent maintenance cost of the equipment is low, and the equipment also retains a good expandability. The equipment has a built-in display screen and a mobile phone display of the excitation gun head, which can operate the software directly on the excitation gun, and visually display the analysis results during the excitation process.

12.The First New Design

It is light in weight and small in size, easy to carry, analysis and measurement is not limited by the site, and can be installed in a matching suitcase and can be brought to the site for material testing anytime, anywhere. At the same time, it can also be easily brought on a high platform to perform inspection work at high altitude and complete outdoor operations.

Operating Requirements

1 Environmental Requirements

The instrument must be placed in an area free from harmful, flammable and corrosive gases. Note: It is forbidden to place this instrument in a chemical analysis laboratory.

Storage temperature: 0 c to 45

2 Power Requirements

- 1)Power supply: 24V/10A.
- 2)Charging power supply: single-phase 220±20V power supply.

3Argon Requirement

1)The purity of argon gas is greater than or equal to 99.999%, the oxygen content is less than or equal to 2ppm, and the content of H₂O is less than or equal to 5ppm. If high-purity argon is not available, an argon purifier is recommended.

2)Argon flow rate: maintain the flow rate of about (0.2-0.3) L/min, and stimulate the flow rate of about 3L/min.

3)Argon control pressure: 0.4Mpa.

4Argon Connecting Pipe

Must use the special argon tube connection provided.

5. Sample Preparation

An angle grinder is used to grind the surface of the tested sample.

The sample to be analyzed must be uniform, free of pores and casting defects, and the preparation must achieve a smooth surface, no oxide layer, no oil stains, and no burrs on the edge of the sample.

6. Standard Sample Requirements

A spectral calibration standard sample will be randomly attached when it leaves the factory to correct the overall spectral drift of the instrument. In addition, users need to prepare standard samples or internal control samples suitable for their own product types for calibration of instrument analysis curves.

Installation and Acceptance

1. Arrival Inspection and Storage

After the buyer receives the goods, they should check whether the outer packaging is in good condition. If the outer package is damaged or the direct reading spectrometer is damaged, the user should keep relevant evidence and notify the seller in time, and cooperate with the seller to make a claim to the responsible party. Otherwise, the buyer is responsible for damage to the goods.

After the instrument is delivered to the buyer, the buyer should store the instrument in an indoor environment with a temperature of 0°C to 45°C and a humidity of 20% to 80%. Before the installation engineer arrives, the buyer cannot open the box without authorization, otherwise, the goods will be deemed as qualified.

2. Installation and Commissioning

Within one month of the arrival of the instrument, the buyer should meet the installation conditions and notify the seller to install it in time.

After the installation engineer arrives, the packing list will be inspected and received by both the packing list in duplicate according to the inspection results and sign them. If there are any missing parts, the installation engineer shall indicate them on the installation acceptance sheet, and immediately notify the seller to make up the mail.

If the buyer still does not have the installation conditions within two months of the arrival of the instrument, or if the instrument is damaged during storage, the seller has the right to ask the buyer to pay additional installation costs, and all costs caused by the damage to the instrument during storage shall be borne by the buyer.

During the installation and commissioning of the instrument, the Buyer shall provide the

installation engineer with convenient working conditions. If the installation engineer can't work and stay at the site for more than 2 days (inclusive) because the buyer doesn't have the installation conditions, the installation engineer can directly return, and the buyer shall bear all the expenses of the installation engineer during his stay, including round-trip transportation, accommodation, subsidies, etc.

Installation Step:

1) Engineers go to the site to install and debug according to the instrument installation procedure.

2) Selection of test standard sample:

According to the content and grade of the samples analyzed by the buyer, the buyer and the seller shall negotiate to determine a set of standard samples that meet the daily analysis requirements of the user (the set of standard samples shall be provided by the buyer or entrusted by the seller to purchase through negotiation) as the standard sample for instrument acceptance. The principle of standard sample selection should be based on the principle of meeting the content range required by the buyer's daily production.

After-Sales and Others

1. Random Technical Data and Software

User Manual

Daily Maintenance Manual

2. Training

During the on-site installation and commissioning of the installation engineer, the buyer shall organize relevant users and maintenance personnel to receive on-site training, and the training time generally does not exceed 3 days. Through the training, the relevant personnel can have a deeper understanding of the relevant technologies of spectral analysis and be proficient in the operation and maintenance of the instrument.

3. After-Sales Service

1. Calculation of Warranty Period

1) The warranty period of the instrument is subject to the contract

2) Unless otherwise specified in the contract, the warranty shall be implemented in accordance with the following periods:

The warranty period is 12 months after the acceptance of the goods or 12 months after the buyer receives the goods, whichever comes first.

The warranty period of general wearing parts is three months from the date of shipment;

The repair warranty period is three months;

Consumables (standard samples, electrode brushes, tungsten electrodes, etc.) are not guaranteed.

3) If there is a quality problem with the instrument that affects its use, the warranty period will be extended accordingly.

2. During the warranty period, the failure under the normal use of the instrument will be repaired free of charge; However, due to the damage to the instrument hardware caused by the operation or improper use of the buyer's operation or maintenance personnel, as well as force majeure factors (war, earthquake, lightning strike, flood, etc.), the seller will charge the cost fee and maintenance fee.

3. After the instrument is reported for repair, the seller shall first confirm the fault phenomenon by telephone and solve the problem. If it is necessary to solve the problem on site, it can rush to the site within 48 hours.

4. Out of Warranty Maintenance:

After the warranty period expires, the seller will continue to provide maintenance service and use consultation. If it is necessary to provide on-site service to users, the seller will properly charge the spare parts fee, service fee and other fees.

5. Supply of Spare Parts and Consumables:

The seller will supply spare parts and consumables to the buyer at preferential prices.

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies




UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Partical
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

