

RamanMicro-3081

Raman Microscope(Mapping)



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

RamanMicro-3081 combines the advantages of the microscope and the Raman spectrometer. Unique conjugate focusing system makes accurate focusing possible.

The RamanMicro-3081 is equipped with a special lens designed for the Raman system, which makes the laser light spot close to the diffraction limit, and the focus information is displayed accurately and directly on the computer through the 5 million cameras. Overcomes the problem of collecting Raman's signal in the ordinary Raman system slightly higher or lower than the actual optimal focal plane, thus improving the spectral quality. RamanMicro-3081 unique software optical path switching part perfectly solves the loss of camera imaging time and realizes the separation of camera imaging and Raman signal collection, thus obtaining the optimal strength.

Meanwhile RamanMicro-3081 used specifically for high performance Raman micro Raman system optimization, sensitivity, signal-to-noise ratio, stability, etc., are the industry leading level, provide strong guarantee for Raman

►► Features:

- True focus, guaranteed more accurate Raman image.
- Ultra-high spatial resolution.
- Unique software control switching path Super high stability.
- Quickly locate and quickly locate the focus position.
- High quality objective lens, spot micron scale.
- 5 million cameras, accurate image.
- Laser wavelength can be selected from 532, 785 and 1064 nm.
- Equipped with a high performance spectrometer.
- The usb 2.0/LAN interface is directly connected to the computer.

►► Application

- Biological sciences
- Forensic analysis
- Material science
- Medical immunoassay
- Agricultural and food safety
- Water pollution analysis
- Gemstone
- Environmental science

RamanMicro-3081 Automated focusing Raman Microscope	
Spectral resolution	4 cm ⁻¹
Spectral range	200-3800、250-2700 cm (Other wave ranges can be customized to a minimum of 50cm-1)
Spectral stability	$\sigma/\mu < 0.5\%$ (COT 8 hours)
Temperature stability	Spectral shift ≤ 1 cm ⁻¹ (10-40°C)
SNR	>5000:1
Detector	TEC refrigeration 2048*64 pixels CCD
SN ratio	>5000:1
Wavelength range	200nm-1100nm
Pixel size	14 μ m*14 μ m
Detector dynamic range	13000:1
Laser center wavelength	785nm (+/-0.5nm)
Microscope camera	3-megapixel /5-megapixel camera
Focusing Mode	Conjugate focusing(True confocal)
Laser power	>550mW (software adjustable)
laser beam diameter	>1 μ m
Laser stability	$\sigma/\mu < \pm 0.2\%$
Laser linewidth	0.08 nm
Communication mode	USB2.0



Fig. 1 RamanMicro-3081 Structure indication diagram

►► 2. Optical performance

►► 2.1 Spectral exhibition




Fig. 2 RamanMicro-3081 collect acetonitrile spectra







Fig.2 RamanMicro-3081 Sers experiment 1 (Left picture is sample, and right picture is Sers Raman spectra)



RamanMicro-3081 Sers experiment 2 (left picture is sample, right picture is Sers Raman spectra)




Fig 3 RamanMicro-3081 Measure Si Raman spectra (500mW, integration time: 1S)




Fig 4 RamanMicro-3081 measure alcohol spectra (500mW, integration time:1S)




Fig 5 RamanMicro-3081 measure diamond Raman spectra (30mW, integration time: 1S)




Fig 6 RamanMicro-3081 measure boron carbide (PN) spectra (200mW, integration time: 2S)

►► 2.2 Raman resolution

►► 2.2.1 Tylenol Raman spectra

Power output: 200 mW
Integration time: 10 s
Filter level 1

Raman spectra of Tylenol showed the resolution condition in the long wavelength region.




Fig 7 Tylenol spectra shows clear 1610/1615 cm^{-1} vibration peak

Power output: 200 mW
Integration time: 10 s
Filter level 1

Raman spectra of Tylenol showed the resolution condition in the long wavelength region.




Fig 8 93# petrol Raman spectra, 723/732/742 cm^{-1} spectral shift is clearly recognized

►► 3Reliability

Fig3.1, Fig3.2 temperature stability is measured by RamanMicro-3081, keep stable above an hour for each temperature node ranging between 5-40°C. Sample measured is acetonitrile, wavenumbers shift $\leq 1\text{cm}^{-1}$ (Fig 3.1), peak top intensity change $< 10\%$ (Fig 3.2)




Fig. 3.1 Wavenumber shift results testing from 5°C to 40°C of five RamanMicro-3081 portable Raman spectrometers




Fig. 3.2 Intensity Variation testing from 5°C to 40°C of five RamanMicro-3081 portable Raman spectrometers

►► Order guide

PN	Wavelength /nm	Power /mW	Wavenumber range/cm ⁻¹	Resolution/cm ⁻¹
RamanMicro-3081-473	473	100	150-4000	7
RamanMicro-3081-532	532	100	150-4000	7
RamanMicro-3081-785-27	785	600	250-2700	5
RamanMicro-3081-785-40			150-4000	6
RamanMicro-3081-830	830	600	150-4000	7
RamanMicro-3081-1064	1064	600	150-4000	10
Available in custom wavelength				

►► Details



Fig 9 branded high stable microscope platform; X, Y, Z-axis precision adjustable; Adjustable knob work smooth, weight up to 5.6 Kg, very stable.



Fig 10 Raman signal high transmission objective, confocal distance length up to 8nm



Fig 11 Power button, button on/off as many as 1,000,000 times, high strength laser cable, signal cable is very strong, and laser indicator can intuitively display operating status.



Fig 12 Simple interface: Raman microscope: power socket +USB2.0 connector

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Particle
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA




NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide