

SIZER-3002

Wet and Dry Laser Particle Size Analyzer



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company


www.analyticalgroup.net

Sizer- 3002 is intelligent full automatic laser particle size analyzer, adopt the most advanced International MIE scatter-ing principle and Fraunhofer diffraction principle, which is integrated model of wet and dry dispersion system,High sen-sitive and High-resolution photoelectric probe system and imported Canon lens optical path ensure the good accuracy and repeatability of particle size distribution,It's widely used in industrial production quality control departments and re-search institutions

Advantages:


Automatic Optical path alignment System:

which is composed of precise four phase hybrid stepper motor, Its inching precision is reach to micron level, make optimum optical paths to ensure accurate and stable test.



Automatic sample dispersion system

The modular design of dry and wet dispersion shortens the circulation pipeline and prevents the precipitation of large particles. It is not only conducive to the uniform dispersion of samples, but also reduces maintenance costs and is conducive to the cleaning and maintenance of equipment.




Wet and dry, unique separated dispersion module

The instrument adopts a modular design for automatic dry and wet dispersion. Customers can freely assemble it according to the test requirements of sample characteristics. It is simple, convenient and fast.




Instrument Software:

The unique unconstrained free fitting technology collects scattering data during the measurement process and has the ability to automatically analyze multi-modal particle groups. It is not constrained by any function. Particles below 2µm are tested in 12 particle size grades with extremely high physical resolution



Test principle:



Application:

Instrument Software:

The unique unconstrained free fitting technology collects scattering data during the measurement process and has the ability to automatically analyze multi-modal particle groups. It is not constrained by any function. Particles below 2 μ m are tested in 12 particle size grades with extremely high physical resolution




Long service time


Laser particle size analyzer, as an analysis instrument, it doesn't have consumable parts except for the stirring parts, it has no transmission parts and no wearing parts; high performance laser, with a long service time of more than 25000 hours, high sensitive photodetectors is a core part, it will not be easily damaged if operated normally; the photodetector array is a key part, as long as it is used properly, it will not be automatically damaged. Therefore, users do not have to worry about the service time of winner instruments at all.




1. Unconstrained free fitting technology can truly reflect the particle distribution



2. Accurate and convenient automatic alignment function.




3. Different refractive index models can be established to make the measurement results more accurate and reliable



Model Name	Refractive Index	Extinction Coefficient
Model 1	1.5	0.0001
Model 2	1.55	0.0002
Model 3	1.6	0.0003
Model 4	1.65	0.0004
Model 5	1.7	0.0005
Model 6	1.75	0.0006
Model 7	1.8	0.0007
Model 8	1.85	0.0008
Model 9	1.9	0.0009
Model 10	1.95	0.001
Model 11	2.0	0.0011
Model 12	2.05	0.0012
Model 13	2.1	0.0013
Model 14	2.15	0.0014
Model 15	2.2	0.0015
Model 16	2.25	0.0016
Model 17	2.3	0.0017
Model 18	2.35	0.0018
Model 19	2.4	0.0019
Model 20	2.45	0.002
Model 21	2.5	0.0021
Model 22	2.55	0.0022
Model 23	2.6	0.0023
Model 24	2.65	0.0024
Model 25	2.7	0.0025
Model 26	2.75	0.0026
Model 27	2.8	0.0027
Model 28	2.85	0.0028
Model 29	2.9	0.0029
Model 30	2.95	0.003
Model 31	3.0	0.0031
Model 32	3.05	0.0032
Model 33	3.1	0.0033
Model 34	3.15	0.0034
Model 35	3.2	0.0035
Model 36	3.25	0.0036
Model 37	3.3	0.0037
Model 38	3.35	0.0038
Model 39	3.4	0.0039
Model 40	3.45	0.004
Model 41	3.5	0.0041
Model 42	3.55	0.0042
Model 43	3.6	0.0043
Model 44	3.65	0.0044
Model 45	3.7	0.0045
Model 46	3.75	0.0046
Model 47	3.8	0.0047
Model 48	3.85	0.0048
Model 49	3.9	0.0049
Model 50	3.95	0.005
Model 51	4.0	0.0051
Model 52	4.05	0.0052
Model 53	4.1	0.0053
Model 54	4.15	0.0054
Model 55	4.2	0.0055
Model 56	4.25	0.0056
Model 57	4.3	0.0057
Model 58	4.35	0.0058
Model 59	4.4	0.0059
Model 60	4.45	0.006
Model 61	4.5	0.0061
Model 62	4.55	0.0062
Model 63	4.6	0.0063
Model 64	4.65	0.0064
Model 65	4.7	0.0065
Model 66	4.75	0.0066
Model 67	4.8	0.0067
Model 68	4.85	0.0068
Model 69	4.9	0.0069
Model 70	4.95	0.007
Model 71	5.0	0.0071
Model 72	5.05	0.0072
Model 73	5.1	0.0073
Model 74	5.15	0.0074
Model 75	5.2	0.0075
Model 76	5.25	0.0076
Model 77	5.3	0.0077
Model 78	5.35	0.0078
Model 79	5.4	0.0079
Model 80	5.45	0.008
Model 81	5.5	0.0081
Model 82	5.55	0.0082
Model 83	5.6	0.0083
Model 84	5.65	0.0084
Model 85	5.7	0.0085
Model 86	5.75	0.0086
Model 87	5.8	0.0087
Model 88	5.85	0.0088
Model 89	5.9	0.0089
Model 90	5.95	0.009
Model 91	6.0	0.0091
Model 92	6.05	0.0092
Model 93	6.1	0.0093
Model 94	6.15	0.0094
Model 95	6.2	0.0095
Model 96	6.25	0.0096
Model 97	6.3	0.0097
Model 98	6.35	0.0098
Model 99	6.4	0.0099
Model 100	6.45	0.01


4. Automatically memorize the last sample test information and display the current test process in real time. Freely customize the display mode and switch between energy spectrum and data display.




5. User-defined analysis parameters, calculating the percentage according to the particle size, calculating the particle size according to the percentage, or calculating the percentage according to the particle size range, so as to meet the characterization methods of particle size testing in different industries



6. Statistical comparative analysis can be carried out for multiple test results, and the difference between different batches of samples, samples before and after processing, and test results at different times can be clearly compared, which has strong practical significance for the quality control of industrial raw materials.



7. Multiple distribution modes: free distribution (closer to the real data of the sample), Rosin-Ramier distribution, logarithmic normal distribution and original data conversion mode (according to the real and accurate measurement of abrasive and flake particles)



8. English language interfaces are supported, and other language interfaces can also be embedded according to user requirements. Multiple formats can be set for file printing and exporting, and BMP image files, Txt documents, Word documents, and Excel documents.

Hardware Future:

1. Dual laser orthogonal beam patented technology

2. All built-in dispersion systems

It avoids the problem of test data distortion caused by long optical path, uneven dispersion, and large particles settling in the pipeline caused by the external dispersion system

3. Omnidirectional Scattered Light Detection Technology


Adding multiple auxiliary integrated photodetectors can effectively collect scattered light at various angles corresponding to the test range, and achieve test accuracy and reliability in the full range

4. Spectrum Amplification Technology

The sensitivity of the probe to the signal is enhanced, and the test range is greatly improved.

5. Converging light Fourier transform patented technology

The large-angle scattered light is not limited by the aperture of the Fourier lens. The optical path is shortened to the shortest, effectively improving the resolution of the instrument; The optical path design principle belongs to the international leading technology.



6. MIE scattering theory

The full range adopts the most advanced MIE scattering theory.



7. Fully automatic alignment system

The precision four-phase hybrid stepping motor is used to auto-matically adjust the optical path and calibrate the optical path at any time, eliminating the deviation caused by manual alignment, and improving the accuracy and stability of the test results from an optical point of view.

Parameter :

Measurement Range		Wet: 0.01-2000µm Dry: 0.1-2000µm.
Sample Injector		Wet Wirmodule-201 Dry: Winmodule-301
Number of Channels		106пов
Accuracy Error		Wet <0.5%, Dry: <1.0% (CRM D50)
Repeatability Error		Wet: <0.5%, Dry: <1.0% (CRM D50)
Light Source		Main laser: High-performance laser X-639nm. P>2mW Auxillary laser. Blue laser A405nm. P>2mW Service life: >50000h
Optical Bench Alignment System		Automatic alignment
Operation Mode		SOP Fully automated operation mode
Dispersion	Ultrasonic	Frequency: f=40KHz, Power p=50W (Time adjustable)
	Str	Revolution speed: 0-3000rpm (Adjustable)
	Circulate	Rated flow: 1L-BL/min; Rated power: 120W
	Dry Method	Turbulent dispersion, automatic sample introduction, one-click testing
	Wet Method	Ultrasonic, mechanical stiming, and built-in circulating dispersion
Sample Tank Volume		600ml
Test Speed		Fastest 10 seconds
Retractive index		A wide range of refractive indices are available,any refractive index can be entered.
Sampling Frequency		Fastest data acquisition rate of 10kHz
Retractive index		A wide range of refractive indices are available,any refractive index can be entered.
Sampling Frequency		Fastest data acquisition rate of 10kHz
Dispersion	Analysis Pattern	Free distribution,R-R distribution,log-normal distribution,and statistical models by category
	Statistical Comparison	Statistical comparison of multiple test results. Users can compare and analyze feet data from different batches of samples, before and after sample processing or at different processing times,effectively guiding product quality and cost control
	Custom Analysis	The system can customize the output of any characteristic particie size between Di and D100,the cumulative percentage of particles larger or smaller than a certain particle size,and the cumulative percentage of a certain particle size range,so as to meet the different requirements of different ndustites for particle size testing and characterization methods
	Output	Typical particle size values, interval particle size distribution (Le differential distribution) data and graphe(histograms or curves)cumulative particle size distribution(le.,integral distribution) data and curves, etc.
	Test Report	Test reports can be exported in various formats such as Word, Excel, Image (bmp),text (txt), and PDF,meeting different requirements for viewing test reports,data reprocessing,and article citation in any situation
	One-click operation	Wet method:Activate the automatic operation mode and complete all operations such as water intake,air bubble removal,dispersion testing cleaning,and data processing with a single button. Dry method: Activate the sutomatic operation mode to complete all operations such as dust collection, air supply, feeding, testing, and data processing with a single click
	GMP	GMP (Optional)
Outer Dimension		Main unit: 870mm*340mm*380mm Wet module: 415mm*280mm*335mm Dry module: 260mm*245mm*296mm
Net weight		Main unit: 32kg. Wet module 15kg; Dry module: 7kg

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



LCMS



Automated Prep-Flash Chromatography system



Maldi TOF



Optima Gas Chromatograph



Flash Chromatograph



DAC Column



GCMS 3068




UHPLC



HPTLC



Ion Chromatograph



Production HPLC



Helium Mass Spectrometer Leak Detector



Column



DLS



Water purification system

▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility

ANALYTICAL FOUNDATION is a Nonprofit Organization (NGO) for the purpose of:



Analytical Foundation


1. Research & Innovation **Scientist's awards & International Conference visit sponsorship / QC Professional Awards:** Quality life is possible by innovation only and the innovation is possible by research only, hence **ANALYTICAL FOUNDATION** is committed to identify such personallities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research /testing/ publication at info@analyticalfoundation.org

2. Improving quality of life by offering **YOGA Training courses, Work shops / Seminars** etc.

3. **ANALYTICAL FOUNDATION** aims to DETOXIFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Naturecure, Health Care, Ara, Events, Camps etc.

▶▶ A portion of profit of **ANALYTICAL GROUP** Company's goes into **ANALYTICAL FOUNDATION** for the aboce noble cause.

▶▶▶ Reach us @



Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice

HPLC Solutions	MultipleLabs	Infinite Multiplelabs LLP	Analytical Bio-Med	Analytical Foundation (Trust)
----------------	--------------	---------------------------	--------------------	-------------------------------

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E :- info@analyticalgroup.net
info@hplctechnologies.com
info@multiplelabs.com

W:- www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide
USA| UK| Australia| Dubai| Malaysia