

STA FTIR- 3403

FTIR-Synchronous Thermal Analyzer



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

The working principle of simultaneous thermal analysis instrument (TGA-DSC) is based on programmed temperature control to measure the relationship between mass of substances and temperature or time. Its main components include heating furnace, balance, temperature controller and recording system.

Operational Principle:

Based on the temperature change of the sample and the process of heat conduction, absorption, radiation, etc. It controls the temperature of the furnace to heat, cool and maintain the temperature of the sample, and then detects and records the temperature change of the sample and the corresponding physical and chemical properties.

Heating process: The sample is placed in a crucible within the heating furnace and is heated or cooled at a predetermined rate through a programmed temperature control system. During this process, the sample undergoes a series of physical and chemical changes such as moisture evaporation, organic decomposition, phase transitions, etc., which can lead to a decrease in the sample's mass.


Quality monitoring: The quality change of the sample is measured and recorded by the high precision balance in real time. The balance collects and transmits the quality change to the host computer in real time, and finally displays the relationship between the quality and temperature or time in the form of curve.

Data recording and analysis: The thermal stability, phase transition temperature, decomposition temperature and other characteristics of the sample can be analyzed by the quality change curve recorded by the recording system.

Application area:

The simultaneous thermal analysis instrument (TG-DSC) can detect both tg/dsc information at the same time, and can measure a wide range of information, including: thermal stability, oxidation reduction, decomposition behavior, corrosion study, decomposition kinetics analysis, melting/crystallization, solid phase transition, crystallinity, glass transition, etc.

Thermal conductivity analysis instrument (TG-DSC) is widely used in the research and development, process optimization and quality control of plastics, rubber, coatings, pharmaceuticals, catalysts, inorganic materials, metal materials and composites.



Plastic materials



metal material



medicine



paint

Structural advantages:

1. Furnace heating: the double-row winding structure of precious metal platinum-rhodium alloy wire is adopted to effectively reduce electromagnetic interference and have better high temperature resistance.
2. Tray sensor: precision processing of precious metal alloy wire, with high temperature resistance, oxidation resistance and corrosion resistance.
3. Power supply and heat dissipation: The power supply system and circulating heat dissipation module are set independently of the host, effectively isolating heat and vibration interference, and ensuring the measurement accuracy of the microthermometer.
4. Heat insulation of the host: The heat insulation design of the host effectively isolates the conduction effect of the high temperature of the furnace body on the chassis and the micro thermometer.
5. Furnace function: Double insulation structure is adopted to achieve higher linearity of temperature distribution; automatic lifting function is provided to quickly cool down and improve operation safety; exhaust gas outlet is set up to support the connection with infrared spectrometer and other equipment to realize online gas analysis.

Advantages of controller and software:

1. The imported ARM processor is used to achieve faster sampling speed and processing speed.
2. The four-channel sampling AD collects the TG signal and temperature T signal.
3. Heating control, using PID algorithm for accurate control. It can be heated in multiple stages and kept at constant temperature (calcium oxalate repeatability test)
4. The software and the instrument use two-way USB communication to realize remote operation completely, and the instrument parameters can be set and the operation stopped through computer software.
5. 7-inch full-color 24-bit touch screen for a better human-machine interface. All calibration of TG can be achieved on the touch screen





(FTIR-TGA)

Precautions:

temperature range	Room temperature-1550°C
temperature resolution	0.001°C
temperature fluctuation	±0.01°C
Heating/cooling rates	0.1~ 100°C/min°C-0.1~-40°C/min
Temperature control method	PID control, heating, cooling, constant temperature
Weighing system measurement range	0.001mg~3g, can be expanded to 10g
accuracy	1µg
sensitivity	0.1µg
DSC range	±1000mw
DSC sensitivity	0.1µw
Atmospheric control	Equipped with three gas mass flowmeters,including two for sample gas switching and one for furnace purge gas, with program control
software	The intelligent software can automatically record the TG curve for data processing, and can switch between TG/DTG/DDTG, quality, and percentage coordinates at will. The software has an automatic adjustment function, which can automatically extend and zoom according to the spectrum display.
Constant temperature time	Arbitrary setting; standard configuration ≤600min
The operating software has the corresponding copyright certificate, and the data test frequency can be selected as real time, 2S,5S, 10S, etc	
The furnace is equipped with automatic and manual lifting modes, which can cool down quickly; 1000°C to 50°C in 20 minutes or less	
External water cooling device, isolation of heat drift effect on weighing system; temperature range-10~60°C	

Test graph:



HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



LCMS



Automated Prep-Flash
Chromatography system




Maldi TOF



Optima Gas
Chromatograph



Flash
Chromatograph



DAC
Column



GCMS
3068




UHPLC



HPTLC



Ion Chromatograph



Production
HPLC



HPLC



Column



DLS



Water purification
system

▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility

ANALYTICAL FOUNDATION is a Nonprofit Organization (NGO) for the purpose of:



Analytical Foundation

1. Research & Innovation **Scientist's awards & International Conference visit sponsorship / QC Professional Awards:** Quality life is possible by innovation only and the innovation is possible by research only, hence **ANALYTICAL FOUNDATION** is committed to identify such personallities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research /testing/ publication at info@analyticalfoundation.org

2. Improving quality of life by offering **YOGA Training courses, Work shops / Seminars** etc.

3. **ANALYTICAL FOUNDATION** aims to DETOXIFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Naturecure, Health Care, Ara, Events, Camps etc.

▶▶ A portion of profit of **ANALYTICAL GROUP** Company's goes into **ANALYTICAL FOUNDATION** for the aboce noble cause.

▶▶▶ Reach us @



 **Analytical**®
Technologies Limited

Subsidiaries



Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice

HPLC Solutions

MultipleLabs

Infinite Multiplelabs LLP

Analytical Bio-Med

Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E :- info@analyticalgroup.net
info@hplctechnologies.com
info@multiplelabs.com

W:- www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide
USA| UK| Australia| Dubai| Malaysia