

Transmission Electron Microscopy Miniature



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Purity analysis from development to commercial process

Transmission electron microscopy (TEM) can provide unique insights when characterizing viral gene delivery platforms such as those based on adeno-associated viruses (AAVs). Morphological characterization can support process development by confirming purity and status of the viral capsids.




►► Example of a downstream process for AAV

- As more and more gene therapy projects move into later clinical phases, manufacturers face new challenges. When developing the commercial manufacturing process scaling up may have impact on the purity profile at different steps. Access to reliable analytical data to base process change decisions on is key to success for a smooth process development.

Impurities or morphology features in AAV samples that can be detected by TEM:

- Broken particles
- Proteasomes
- Protein or membrane-based cell debris
- Residual DNA
- Capsid content (full/empty)
- Aggregates
- Remnants of helper virus



►► Miniature in process development


- Experience from providing an EM service for gene therapy vectors, and of the type of features that need to be characterized frequently and in high volume, has inspired AS to develop Miniature. Miniature is a low voltage system that can be placed in any standard lab, close to your process. Miniature provides bioprocess workers with their own in-house solution that enables non-experts in electron microscopy to generate quantitative data and images in a few hours - data that can support decision making in process development.

Miniature:


- Automates imaging, particle detection and classification
- Visualizes your product and confirms critical quality attributes
- Turns visual evidence into quantitative data
- Generates instant results that allow you to adjust the process within hours

►► Monitoring each step in the purification of AAV


Negative staining (TEM) is a powerful method for assessing process- and product-related impurities. It is a quick and easy preparation technique, supported by Miniature, in which the specimen is embedded in a layer of heavy metal salt solution. The use of heavy metal stain significantly increases the contrast in the images.



Sample after 1st
purification step




Sample after 2nd
purification step




Sample after final
purification step


►► The types of contaminants that may appear




(A) Large size impurities, 4 μm



(B) Small size impurities, 2–15 nm.



(C) Aggregates



(D) Residual DNA impurities, 2 nm thickness


Undesired sample characteristics to consider when optimizing a process

- Aggregation
- Failure to remove host-cell debris or residual DNA
- Loss of particle morphology or integrity


All these types of contaminants can easily be identified and quantified by using the automated image analysis software tools in Miniature.

►► Automated analysis of purity and level of debris


Miniature automatically images, detects and classifies particles and debris. Since automation enables detection and analysis of a far larger number of particles than is feasible with manual operation, the results are more consistent, statistically significant and reproducible.




Undetected particles



Detected primary particles and debris




Particles detected, classified and quantified using Miniature.



Manual imaging on conventional EM was not performed in its entirety, for practical reasons. In order to compare the time required to obtain data of similar accuracy using manual handling, the time for data analysis of 20 images was extrapolated to estimate the time that would be needed for 504 images.


►► Proteasomes could be hiding in your AAV sample

Comparing Miniature analysis with DLS size distribution analysis on the same sample demonstrates that the quantitative data based on image analysis by Miniature can reveal the presence of proteasomes where DLS gives only one broad peak.



►► Size distribution analysis to reveal aggregates and debris




Particles present in the sample that deviate from the expected size and shape can be the undesired outcome of an un-optimized process step. In this example a large aggregate as well as smaller size debris are detected and automatically measured.




Outliers represent potential aggregates. By clicking on the outlier Miniature takes you to the image and you can determine whether it's a true aggregate or not.

►► AAV capsid integrity variations observed with Miniature

Negative-stained samples of AAV particles display variations in staining patterns that provide a clear indication of capsid integrity status. Intact AAV particles appear as lacking internal staining as opposed to broken AAV particles that exhibit internal staining because of presumed disrupted virus capsid structure. Particle doublets and triplets that can occur under certain conditions can also be detected.



►► Using Miniature to compare the capsid status in two different samples




In this Miniature study, the integrity of two different samples of purified AAV particles was automatically quantified and compared.

The study showed that AAV Sample 1 had a larger portion of intact particles (A) compared with Sample 2 (99% intact particles compared with 97% in Sample 2. (B) No doublets or triplets were detected in Sample 1, whereas 42 doublets or triplets (0.125%) were detected in Sample 2.

►► Linearity proves accuracy of TEM full/empty analysis

TEM is a method by which the biological specimens are vitrified by rapid freezing in liquid ethane, embedding the molecules in a layer of amorphous ice. TEM-method does not involve using stains or other chemicals, so therefore samples prepared in this way preserve their native structure. The particles observed with TEM-method with internal dark staining represent particles containing DNA. TEM-method in combination with AS Analyzer Software (AS) has been validated for quantification of percentage full AAV analysis. In this example three different samples are analyzed: empty (S1), filled (S3) and a third sample that is a mix 1:1 of sample S1 and S3 (S2). The TEM-method shows the expected linear relationship between particle content and staining.




The quantitative image analysis of 3 imaged samples defined as empty (S1), 51% filled (S3) and a mix of S1 and S3 (S2) performed with AS is shown in the graph in figure (A). The imaging results are shown in figure (B) and (C) below.

►► Validated AAV full/empty analysis service

Virus particle packaging analysis with TEM-method in combination with AS is recommended for batch release testing and QC of AAV-based gene therapy products. The automation that AS provides allows AS EM services to deliver high quality images and timely and cost-effective accurate quantitative data.

Scatter plot showing 2 distinct populations of particles; full in red and empty in blue. Particles that don't clearly belong to any of these groups and which appear between the red and blue clusters are identified as uncertain.



►► Analyzer software

analyze transmission electron microscopy (TEM) images of nanosized particles.
enables reproducible and semi-automated particle detection and classification. For QC testing of capsid percentage full/empty of AAV samples AS offers services using TEM-method and image analysis that can be validated according to GMP.



Detect

The software works by algorithmically identifying particle characteristics based on selected criteria, such as the size or shape.



Explore

Each detected particle is stored together with the meta data. You can use the collected data to qualify particles into separate classes or groups, such as filled or empty.



Report

After the analysis generates a report of the results with particle data presented in histograms, with statistics, together with representative microscopy images.



Compliance with GMP

The software is 21 CFR part 11 compliant and has access control, electronic signature, full traceability and audit log. suitable for regulatory controlled environments that work under GMP.

Full traceability

- Database architecture for safe archiving
- Instant import of images from live microscope session
- Pre-defined roles for access and right to process data and sign off
- Electronic sign off



HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph




Atomic Absorption
Spectrophotometer



Liquid Partical
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, ASna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide