

Genius SIZER - 3230

Wet and Dry Laser Particle Size Analyzer



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► introduction:

Genius SIZER 3230 intelligent full automatic wet&dry laser particle size analyzer adopt laser diffraction theory(Mie and Fraunhofer diffraction), measure size is from 0.01 μ m to 2000 μ m(dry 0.1 μ m-2000 μ m), Which offer reliable and repeatable particle size analysis for a diverse range of applications.It use dual-beam& multiple spectral detection systems and side light scatter test technology to significantly improve precision and performance of test, It's the prior choice for industrial production quality control departments and research institutions.

Model Name		Genius SIZER -3230
Standard		ISO 13320-1:2009,GB/T19007-2016,Q/0100JWN001-2013 Compliance with 21 CFR Part 11
Principle		Laser diffraction principle
Analysis		Mie and Fraunhofer scattering
Detector Arrangement		Log-spaced array, test angle from 0.015 degree to 145 degree
Measuring Range		Wet:0.01 μ m-2000 μ m Dry: 0.1 μ m-2000 μ m
Channels Number		Wet:127 pcs Dry:100 pcs
Accuracy error		Wet<1% Dry<1% (CRM D50)
Repeatability error		Wet<1% Dry<1% (CRM D50)
Light source		Dual lens, He-Ne laser P>3.0 MW (λ = 632.8nm) Auxiliary green semiconductor laser (λ = 532.8 nm) P> 2.0MW ,
Laser Safety		Class 1
Wet dispersion	Ultrasonic	Frequency:40KHz Power:60W, Time: \geq 1S
	Stir	Revolutions Speed: 0-3000RPM (Adjustable)
	Circulate	Rated Flow:30L/min Rated Power:70W
	Sample tank	Volume:1000mL
	Micro-Sample cuvette	Volume: 10mL (Available)
Dry dispersion		Dry-turbulence dispersion patent, normal shock wave shear technique
Feeding Speed		Adjustable (Variable speed knob)
Operation Mode		Full automatic / manual control, freely choose
Dispersion medium		Compressed Air, pressure: 0 to 6 bar
Optical bench alignment system		Full automatic, precision is up to 0.2 μ m
Test Speed per time		Wet: <2 Min Dry : <1min Typical measuring time<10S
Outer dimension		L104cm×W44cm×H54cm
Net Weight		70Kg

►► Main Features:

1) Wet and dry sample dispersion system Integrated Design

Genius SIZER 3230 intelligent laser particle size analyzer which is 1st set laser particle size analyzer integrated wet and dry dispersion test in one in China, it successfully resolved the problem of dry and wet technology integration, realize one key to switch, apply to test all the particle size distribution from 0.01 μ m-2000 μ m particles.

2) Intelligent full automatic operation system and manual operation mode, freely choose.

With intelligent automatic mode of operation, to achieve a key test, as long as according to the prompt addition of sample, click the "test", all process will be complete automatically, not only reduce the testing workload, but also eliminate the interference of human factors, to further improve the accuracy and authenticity of testing results.

Wet method: A key to complete water-supply, ultrasonic dispersion, circulation, testing, cleaning, data record, data analysis, save and print are automatically completed, only take 2 mins.

Dry method: A key to complete the dust collecting, air supply, feeding, testing, data processing and other operation, take 1 min.

3) Stable and unique optical path system patented technology

Converging light Fourier transform path system, enables scattering light be not restricted to the lens aperture limit, and Dual-laser orthogonal light technology make use of the semiconductor auxiliary laser extend the test angle from 0.015 degree to 145 degree, ensure receive all the angles of signals.

4) Automatic Optical path alignment System,

The precision of four hybrid stepping motor in the automatic system of optical components, micro precision of 0.1 μ m, the instrument of light path is always at its best to eliminate manually on the light path and the troubles and difficulties but also enhance the accuracy and stability of test results.

5) Full built-in Sample dispersion system.

Auto wet dispersion system, SOP realize one key operation. set mechanical stirring, ultrasonic dispersion, and circulation path in one, It ensures particles uniform dispersion and distribution, avoids many bad phenomenon, such as uneven distribution of particles, large particles deposit because of the long outer dispersing system tube, And it guarantees the representativeness of test result. For dry dispersion system, Turbulence dispersion patented technology and Normal shock shearing effect, make particles sufficient dispersion, ensure good test.

►► Software Function::

1, Analysis Mode

Free Distribution, R-R Distribution, Logarithm Normal Distribution, Mesh number classification etc. meet different demands of particle size statistic in different industries.

2,Statistic Method

Volume Distribution, Quantity Distribution

3,Statistic Comparison

Statistic Several Testing Results to compare and analyze Get difference by compare test result of different batches of samples, samples before and after processing, and different time.

Have great practical significance to industrial raw materials quality control

4,User-defined Analysis

Figure out percentage according to the particle size

Figure out particle size according to the percentage

Figure out percentage according to the particle size range

Meet demands of representation of particle test in different industries.

5,Test Report

Word, Excel,Photo(Bmp), Text etc.


6,Multiple language Support

Chinese&English (Others are available)

7,Intelligent Operation Mode

Automatically control water inflow, dispersion,test and analysis.

Better Repeatability after remove human-factor



►► Application:

Genius SIZER 3230 widely used in cement, ceramics, medicines, lotions, paints, dyes, pigments, fillers, chemicals, catalysts, drilling mud, abrasives, lubricants, coal, sediment, dust, cells, bacteria, food additives, pesticides, explosives, graphite, photographic materials, fuel, ink, metal and non-metal powder, calcium carbonate, kaolin, coal slurry and other powdered materials.

►► Adopt Patents Technology:

- Optical bench design is protected by patent No.- ZL 2014 2 0378380.8,
- Three dimensional-optical bench alignment system is protected by patent No.- ZL 2013 2 0835882.4.
- MIE scattering principle application patent No.- ZL 2013 2 0812021.4.
- Dry particle size analyzer full sealed sample cuvette application is protected by patent No.- ZL.2011 2 0267646.8.
- Dual laser beam orthogonal application is protected by patent No.-ZL 2007 2 0025702.0
- Powder dispersion pump design application is protected by patent No.-ZL 2007 2 0018648.7

►► Packing List:

Items	Unit	Quantity	Items	Unit	Quantity
Power cable	pcs	2	Quality card	pcs	1
COM cable	pcs	2	medical silica tube	meter	1.5
Oil-free air compressor	set	1	2A fuse	pcs	5
vaccum cleaner connector	pcs	1	15A fuse	pcs	5
vaccum cleaner	set	1	PTFE tape	box	1
PU tube	meter	5	USB transform cable	pcs	2
filter	pcs	1	Silica tube	meter	1.5
wrench 8/10	pcs	1	Brush for sample tank	pcs	2
wrench 17/19	pcs	1	Quartz sample cuvette	pcs	1
sampling spoon	pcs	2	Sample cuvette brush	pcs	2
Balloon for clean dry sample tank	pcs	1	Tool box	pcs	1
Cleaning paper	pcs	3	standard sample D50= 2.1μm	bottle	1
inner hexagon screwdriver	pcs	1	standard sample D50=110μm	bottle	1
small flathead screwdriver 3*75	pcs	1	Manual	pcs	1
Small flathead screwdriver 6*38	pcs	1	Warranty card	pcs	1
flathead screwdriver 5*75	pcs	1	Wrench 5.5/7	pcs	1
Software CD	pcs	1	Genius SIZER 3230A	pcs	1

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.



 **Analytical**[®]
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide