

# ZSIZER-3300

Zeta Potential & Particle Analyzer



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

**Analytical Technologies Limited**

An ISO 9001 Certified Company

[www.analyticalgroup.net](http://www.analyticalgroup.net)

**Measurement of the charged state on particle-surface changing with solvents**

**Effective evaluation of the dispersion stability of colloidal solution**

**Isoelectric point measurement from pH response**


**Multifaceted evaluation of colloidal particles with zeta potential and particle size measurements**

## ►► **Diverse information and reliable data by microscopic observation**

In ZSIZER, electrophoresing particles are automatically tracked by the combination of microscopic electrophoresis and advanced image processing. The zeta potential values of individual particles are automatically calculated from their mobility.

The ZSIZER measurement process is not so-called "black box" but observable at any time on a monitor.

In addition, ZSIZER-3300 can measure the particle size by automatic tracking of the particle Brownian motion (by optional function)



## ►► **Application Fields**

### **Characterization of particles / Dispersion stability evaluation of colloidal systems**

**Environment** : Water treatment, waterworks, flocculation agents, dispersants, flocculation control, microorganisms, planktons, asbestos, bubbles, soil, beneficiation technology, and etc.

**Nano-materials** : Functional material development, recording materials, pigments, ceramics, catalytic agents, polymers, carbon black, carbon nanotubes, inkjet inks, aqueous/non-aqueous solvent paints, fuel cells, coating materials (for automobile parts / electronic parts), paper manufacturing, surface active agents, and etc.

**Life Science** : Erythrocyte, biological cells, protein, DDS, liposome, vesicle, drug discovery, and etc.

## ►► **Zeta Potential Measurements**

### **Automatic Measurement Modes**


Four different measurements can be selected from automatic measurement modes:

**1. Zeta potential & Histogram** : Standard measurement method: the zeta potential of particles on the stationary level is calculated on the set measurement conditions. The measurement result is displayed in the histogram.


**2. Flow speed distribution in the sample cell** : Flow velocity (FV) at each position in the cell is measured and the FV distribution is graphed. You can check the nonuniformity of electroosmotic flow so that measurement can be performed at any position depending on the shift from the stationary level.

**3. pH responsiveness measurement (Isoelectric point measurement)** : By inputting the pH of dispersion medium, isoelectric point and pH responsiveness can be analyzed from the change in zeta potential.


**4. Settling / Rising velocity measurement** : In settling/rising velocity measurements, the particle tracking direction is defined as the Y-axis (vertical) direction. The settling velocity of aggregates or coarse particles and the rising velocity of bubbles or hollow particles can be measured.




Zeta potential histogram



Flow velocity distribution in the sample cell



pH responsiveness (Isoelectric point)



Settling / Rising velocity measurement

## Real-time measurement / Archive measurement

In the real-time measurement, you can observe and measure the sample in the measurement cell simultaneously. The image files are saved every one second. Based on image archive, archive measurements can be performed in the different measurement conditions.

You can also double-check the measurement result of specific particle by archive measurement which enables you to judge the existence of impurity and provides you enhanced measurement accuracy.

Besides automatic measurement mode, manual measurement mode can be selected. Manual measurement is suitable for :

- Particles with heady flickering due to irregular shape.
- Particle size is widely distributed in the sample so that the data for each particle size need to be acquired.

The temporal change of zeta potential can be analyzed by using the data of individual particles.

## ►► What is zeta potential?


Zeta potential is defined as the magnitude of charge at the 'sliding surface' in the 'ion diffusion layer' around a colloidal particle.

Colloidal particle dispersed in a solution is positively or negatively charged. In order to keep the particles electrically neutral in the solution, ions with opposite polarities build up around the particles and form 'ionic stationary bed layer'.

Outside of it, there is an 'ion diffusion layer' in which ions with same and opposite polarities coexist. Since ions slide from the outside of the ionic stationary bed layer, this boundary is called 'sliding surface'.

When the zeta potential value of positively or negatively charged colloidal particles is large, the dispersion stability of the sample is good. On the other hand, when the zeta potential is close to zero, the particles have tendency to aggregate and the state of the particle system is unstable.

Zeta potential provides you a guideline for the dispersion stability, aggregability, sedimentation property of small particles. Also, the charged state of surface functional group in a dispersion medium can be evaluated



## ►► ZSIZER sample measurement cells



Standard cell



20 mm cell for non-aqueous solution



Temperature controlled cell for aqueous solution

## ►► Particle size distribution measurements by Brownian motion (option)

Brownian motion is a random kinetic phenomenon which is caused by the irregular collision of solvent molecules being in thermal motion in a colloidal solution. For large size of particles, they need more energy for the particle movement. The larger surface area they have, the more molecular collisions occur in various directions. Since the effects of the molecular collisions cancel each other, the Brownian motion becomes small for large size of particles.

In this way, the particle size and the particle movement due to Brownian motion are related with each other. By applying the Einstein-Stokes equation, the particle size can be calculated from the travel distance of a particle.

In ZSIZER-3300, the random movement of particles is tracked by the image analysis. Then the particle size and its distribution are automatically calculated from the travel distance of a particle and various parameters (such as solution temperature and absolute viscosity).

## ►► Specifications (ZSIZER-3300)

Measurement method	Microscopic electrophoresis (zeta potential)
Zeta potential measurement range	-200~200 mV
Mobility	-20~20 cm <sup>2</sup> /sec·V
Measurable particle size (*)	0.02 μm ~ 100 μm (for zeta potential measurement) 0.05 μm ~ 100 μm (for particle size measurement)
Supply voltage	0~350 V DC *Applicable voltage to electrodes (compatible with external power supply)
Light sources	LED / Laser diode
Source position (Irradiation method)	Scattered / transmitted light method
Camera	CCD black-and-white video camera
Objective lens	x10 objective lens (option: additional lens with another magnification)
Video output	Video signal NTSC
Cell stage	Digital display in 0.001 mm increments *Accuracy 0.01 mm
Measurement cell	Standard measurement cell for aqueous solution (option: cell for non-aqueous solution)
Size / Weight	300 (W) × 600 (D) × 332 (H) 25 kg
Power supply	100 V 1 A 50/60 Hz

## HPLC Servicing, Validation, Trainings and Preventive Maintenance :

**HPLC Servicing** :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

**Trainings** :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

**AMC's/CMC** :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

**Validations** :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

**Instruments** :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



## About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

## Our Products & Technologies



UV/VIS  
Spectro 2080+  
Double Beam



Infra FTIR



Optima Gas  
Chromatograph  
3007



Optima Gas  
Chromatograph  
2979 Plus



Flash  
Chromatograph



Atomic Absorption  
Spectrophotometer



Liquid Particle  
Counter



Optical Emission  
Spectrophotometer



DSC/TGA



Semi Auto Bio  
Chemistry Analyzer



HEMA 2062  
Hematology  
Analyzer



Micro Plate  
Reader/Washer



URINOVA 2800  
Urine Analyzer



Total Organic  
Carbon 3800



Fully Automated  
CLIA



NOVA-2100  
Chemistry Analyzer



PCR/Gradient PCR/  
RTPCR



TOC  
Analyzer



Laser Particle  
Size Analyzer



Ion Chromatograph



Water purification  
system



## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



**Analytical**  
**Foundation**

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [Info@analyticalfoundation.org](mailto:Info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## Reach us @



**Analytical**<sup>®</sup>  
Technologies Limited

HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W. [www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)  
[www.ais-india.com](http://www.ais-india.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide