



## Anesthesia Machine AM-3800

### ▶▶▶ Specification:

- A. Pipeline attachments for Oxygen, Nitrous oxide and compressed air.
- B. Yoke assembly for Oxygen and Nitrous oxide with pin index system.
- C. Durable main switch to put the machine in the on or off position.
- D. Cascade double tube bobbin type flow meters for oxygen and nitrous oxide and single for air.
- E. Safety features like:
  - F. Minimum oxygen flow of 50ml/min or more even when the machine is in ON position.
  - G. Provides 25% or more of oxygen when an anesthetic gaseous mixture is in used.
  - H. Provided with mechanical hypoxic guard.
  - I. Oxygen flush with a flow rate of more than 35L/min.
- J. Machine is supplied with one Sevoflurane vaporizer. The Machine is capable of holding two Seletatec vaporizers (Isoflurane, Sevoflurane) simultaneously.
- K. CO<sub>2</sub> absorber system with the following features:
  - L. Single/Double canister
  - M. Autoclavable
  - N. Canister capacity of 1.2kg or more.
  - O. It is possible to bypass the canister if removed during clinical cases to change sodalime.
- P. APL valve assembly and Bag mount is conveniently placed.
- Q. Independent port for Open circuit.
- R. Provided with two or more drawers.
- S. Machine has a good quality handle and castors
- T. Machine will be supplied with all standard accessories.

## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



**Analytical**  
**Foundation**

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [info@analyticalfoundation.org](mailto:info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## Reach us @



 **Analytical**<sup>®</sup>  
Technologies Limited

HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W. [www.ais-india.com](http://www.ais-india.com)  
[www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.