

BOD Incubator 3000 Series



FEATURES

This is special type of incubator designed to meet the requirements for incubation of biochemical Oxygen demand (BOD) on water and sewage. Also recommended for preservation of vaccines at low temperature and supply of endymatic digestion Processes and other tests demanding constant and uniform low temperature. The unit is highly suitable for work at Temperature below as well as above the ambient. Our BOD incubator is being also used as Environmental Study and Growth Chambers by incorporating Lighting & cooling and heating systems. Inside walls are of Stainless Steel And Out Side of mild steel(SS-304).Tripled walled as above finish, gap filled with Glass wool/PUF insulation Illumination : Chamber illuminated with lamp, illumination switch on control panel at extra Cooling : Hermetically Sealed Compressor & System mounted at top/bottom.

SPECIFICATION

- Temperature range: 5-60°C uniform cooling with cooling embedded in the walls.
- Resolution : 0.1°C.
- Temp. Accuracy : 0.5°C.
- Capacity Cft / Ltr : 21 / 600
- Inner Size (WxDXH– CM) : 80 X 75 X 100
- Control : Microprocessor Based PID Controller with software lock, 0.1°C resolutions.
- Size Inside Chamber Trays
- Forced air circulation ensure uniformity of the temperature.
- Three perforated removable Aluminium or stainless steel shelves.
- Suitable to work on Volts, AC 50 Cycles, use of transformer to regulate Line advised.
- Trays : 3
- Stabili.KVA : 5 .

EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specification of the catalogue without prior notice.