

BOD INCUBATOR

3000 Series



EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Specifications

- BOD Incubator Temperature Range: 10.0 C to 60.0 C.
- Cooling Incubator Temperature Range: 2.0°C to 8.0°C.
- Accuracy / Uniformity: $\pm 0.5^{\circ}\text{C}$ / $\pm 1.0^{\circ}\text{C}$, Resolution: 0.1°C .
- Temperature Display: Digital LED 3 ½ digit, Control: Microprocessor based.
- Temperature Sensor: PT – 100.
- Cooling: CFC free Copeland compressor utilizing R 134A refrigerant.
- Air Circulation: Flange motor with impeller / blower, Insulation: PUF.
- Chamber Illumination: Fluorescent light with door switch.

►► Construction

- Double walled, gap between walls filled with PUF insulation.
- Standard model: Stainless steel mirror polished chamber and heavy gauge powder - coated CRC exterior.
- GMP model: Stainless steel mirror polished chamber and exterior of dull finish stainless steel.
- Double door: Inner full viewing door and outer door insulated with lock and key arrangement.
- Unit mounted on castor wheels.
- Control panel with electrical switch gear and refrigeration system mounted on top for easy service accessibility.
- Trays: SS wire mesh heavy duty.

▶▶ **Servicing, Validation, Trainings and Preventive Maintenance :**

- Servicing : We have team of service engineers who can attend to any make of instrument promptly @the most affordable cost.
- Trainings :We also take up preventive maintenance to reduce downtime of instrument's Trainings.
- AMC's/CMC : We offer user training both in-House and at customer sites on instrument principles, operations, troubleshooting.
- Validations : We have protocols for carrying out periodic Validations as per GLP/ GMP/USFDA norms.
- Instruments : We offer instruments/Renting Services Modules like pumps, detector etc. on Rent.



▶▶ **About Analytical Technologies**

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO :9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

▶▶ **Our Products & Technologies**




UV/VIS Spectro 2080+ Double Beam



Infra FTIR



Optima Gas Chromatograph 3007



Optima Gas Chromatograph 2979 Plus



Flash Chromatograph



Atomic Absorption Spectrophotometer



Liquid Partical Counter



Optical Emission Spectrophotometer



DSC/TGA



Semi Auto Bio Chemistry Analyzer



HEMA 2062 Hematology Analyzer



Micro Plate Reader/Washer



URINOVA 2800 Urine Analyzer



Total Organic Carbon 3800



Fully Automated CLIA



NOVA-2100 Chemistry Analyzer



PCR/Gradient PCR/ RTPCR



TOC Analyzer



Laser Particle Size Analyzer



Ion Chromatograph



Water purification system

▶▶ **Regulatory compliances**



▶▶ **Corporate Social Responsibility**

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

▶▶ **Reach us @**



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalbiomed.com

W: www.ais-india.com
www.analycalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents /technical specifications of the catalogue without prior notice.