

# Centrifuge 3000+



## Technical Specification:

- Maximum tube size:- 100ml
- Maximum capacity:- 400 or more
- Maximum speed:- 20000rpm.
- Temperature Range:- -10°C to +40°C
- Speed holding accuracy:- S100RPM
- Display:- LED display (Microprocessor control)
- Rotor:- 8 places of 25 ml each
- Fixed - Angle Rotor :- 11,500g
- Dimensions (HxWxD):- 216x231x292 mm
- Weight:- Approx, 6 kg
- Frequency:- 50-60 Hz
- Induction motor with frequency drive
- Automatic unique safety rotor identification system
- Safety lid interlock to prevent cover opening during centrifugation.

ISO  
9001

CE

## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



# Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [info@analyticalfoundation.org](mailto:info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## Reach us @



HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W: [www.ais-india.com](http://www.ais-india.com)  
[www.analycalgroup.net](http://www.analycalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide

Note: Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.