

Deep Freezer DF-2100



▶▶▶ Features:

- Upright type(Vertical)
- Inside temperature ranging from -10°C to -40°C can be set freely
- Controlling precision is 1
- Stainless steel interior, Outer is painted steel panel
- Front opening lockable door with full height handle
- 110mm foaming insulation 1 set
- Inner 3 shelves Stainless steel with inner door
- Imported Danfoss compressor, 2-4 hours cooling temperature come to -40 steady
- Fan electromotor with internationally famous of Germany EBM
- Microprocessor- based temperature control system, digital temperature display.
- Keyboard lock and password protected configuration page
- Audible and visual alarm
- High or low temperature alarm
- Power failure
- Low battery alarm
- Door ajar Alarm
- Filter blocking system failure alarm with temperature printer, record 7 days data.
- 72 hour battery back up for printer and controller display.
- 4 units casters for easy handling
- Mixture gas with R134a, R404a, R23, R24, N2, R600a
- Certificate ISO 9001, IOS 14001, ISO 13485
- Power Supply 220V/50Hz 1 phase, can change as 220V 60Hz or 110V 60Hz

▶▶▶ Optional:

- Chart recorder, CO2 back up, Storage racks, Boxes

▶▶▶ Specification:

- Capacity: 105 litres
- External Size (WxDxH)MM: 680x580x935
- Inside size (WxDxH) MM: 480x380x665
- Input power: 120W
- Weight: 38kg



▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

▶▶▶ Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.ais-india.com
www.analycalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide