

# 3000 Plus

## Pharma/laboratory Refrigerator



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

**Analytical Technologies Limited**

An ISO 9001 Certified Company

[www.analyticalgroup.net](http://www.analyticalgroup.net)

## Features:

- Laboratory refrigerators are specially designed cold storage cabinets that are used to store variety of vaccines and medicines (human and veterinary) at prescribed low temperature.
- Full View Inner Glass Door, CFC Free Refrigeration.
- Outer body made of GI/MS Powder Coated OR Stainless steel 304 grade & Inner chamber made of Stainless steel 304 grade.
- Laboratory refrigerators are designed with temperature range 2°C to 8°C. Temperature controlled by PID temperature controller with PT100 sensor for temperature settings.
- This controller displays both Set Value (SV) and Process Value (PV).
- Air circulation is mandatory in any laboratory refrigerator to maintain uniform temperature.
- Lab refrigerator is fitted with illumination facility; a LED lamp is fitted inside the chamber, which is door operated. When door opens, it switches ON automatically.
- Hinged type solid door or glass door, and fitted with reliable door handle with spring latch or pull action clamp.
- Fitted with LED lamp inside the chamber.
- Supplied with removable shelves.
- Fitted with castor wheels.
- Power 220-230 Volts AC Supply.

## Inner Chamber Size & Capacity

- (a) 610 x 450 x 410 mm (Capacity: 113)
- (b) 845 x 500 x 400 mm (Capacity: 173)
- (c) 870 x 570 x 560 mm (Capacity: 282)
- (d) 900 x 650 x 570 mm (Capacity: 342)
- (e) 700 x 900 x 650 mm (Capacity: 421)
- (f) 980 x 640 x 700 mm (Capacity: 450)
- (g) 1100 x 700 x 700 mm (Capacity: 600)

## OPTIONAL

- Digital data logger/ USB port
- Circular chart recorder

## HPLC Servicing, Validation, Trainings and Preventive Maintenance :

**HPLC Servicing** :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

**Trainings** :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

**AMC's/CMC** :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

**Validations** :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

**Instruments** :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



## About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

## Our Products & Technologies



UV-TS 2080



Infra FTIR



Optima Gas  
Chromatograph  
3007



Optima Gas  
Chromatograph



Flash  
Chromatograph



Atomic Absorption  
Spectrophotometer



Liquid Particle  
Counter



Optical Emission  
Spectrophotometer



STA



Semi Auto Bio  
Chemistry Analyzer



LCMS



GCMS



XRD



Total Organic  
Carbon 3800



Fully Automated  
CLIA



NOVA-2100  
Chemistry Analyzer



PCR/Gradient PCR/  
RTPCR



TOC  
Analyzer



Laser Particle  
Size Analyzer



Ion Chromatograph



Water purification  
system

## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



**Analytical**  
**Foundation**

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [Info@analyticalfoundation.org](mailto:Info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## Reach us @



Scan Me

HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W. [www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)  
[www.ais-india.com](http://www.ais-india.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.