

Automatic Slide Staining Machines

Automatic Slide Staining Machines not only provides high degree of flexibility and high specimen throughput in Histopathological Laboratories but also maintains maximum **Quality and Safety Standards**

AAS-2700 AUTOMATIC SLIDE STAINING MACHINE (12 STATIONS): A simple, sturdy and economical model of Automatic Slide Staining Machine with 12 stations, recommended for hema-tology and also for staining of blood films and bone marrow spread etc.



▶▶▶ Key Features

Compact Bench Top Model

Perfect H & E Staining for Routine Application

High Specimen Throughput

Minimum Reagent Exposure

Easy to clean

Simple to operate

Highly recommended for Hematology and also for staining of blood Films and Bone marrow spreads etc.

▶▶▶ Technical Specifications

Operation	Automatic
Staining Troughs	12
Trough Capacity	12 Sides
Trough Volume	450 ml
Inlet Outlet Trough	1
Slide Carriers	1
1 Hr Calibrated Timing Discs	6
1/2 Hr Calibrated Timing Descs	x
Notch Cutting Plier	-
Dust Cover	-
Dimensions in mm	560x560x370
Power Supply	110 / 220 V@50/60 Hr
Timing Seqence	12 Stage Programmable
Timing Step	1 Minute



▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

▶▶▶ Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide