

Birthing Table BT-5200



▶▶▶ Specifications:

(A) Description of Function:

1. Delivery bed is used for Baby Delivery and incorporate ideal blend of the patient's individual requirements on comfort and the professional needs of the delivery team, focusing on the aesthetic and functional design of the entire product.

(B) Operational Requirements:

1. Delivery bed is supplied with all accessories.

(C) Delivery Bed should have following essential specifications:

1. It have control devise for making height (44cm to 90cm) and back adjustments.[manual as well as remote control].
2. It have collapsible side rails
3. It have three sectional mattress and seat section should have large perineal cut. The mattress thickness is 50mm or more.
4. Head board and food section can be detached or slides and stores under the bed.
5. Wheels (dia- 6" or 8") provided with locking system.
6. Retractable foot section with indication for locking, so as to convert bed into table.
7. Infusion rods which have adjustable heights, quick release and attaches to all corners of bed.
8. Adjustable leg rests available as an accessory
9. Push grip handles
10. Sliding stainless steel bowl at perineal part of table
11. Catheter bag holder which can be attached on either side of bed
12. Able to give trendelenburg, reverse trendelburg and 60 degree sitting position both mechanically andelectronically.
13. Adjustable foot supports for nursing staff
14. Easy to clean, sterilize (especially blood stains) and maintain
15. Frame is of epoxy powder coated steel
16. Dimensions - Length: Minimum 180 cm and width: Minimum 75 cm
17. Pelvic tilt: - 15 degree.
18. Easy slide calf supports swing into correct positional lock with single lever.
19. CPR release.
20. Weight capacity: 200 (Approx)

▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



**Analytical
Foundation**

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

▶▶▶ Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide