

Blood Weighing Scale Monitor BWS-3000-Series



Special Features:

- LED display, displays the weight and volume with an accuracy of 1 gm/ml.
- Helps better balancing of refrigerated centrifuge.
- Easy conversion of weight to volume.
- Tare provision to account for the weight of the blood bag
- Compact model BWS molded body & Auto calibration

Specification:

Model	BWS-3000-Series
Input voltage	230VAC
Power supply	SMPS
Transducer	Load Cell
Weighing Range	3 Kg
Accuracy	1gm/ml
Calibration	Auto
Display	LCD with back light
Weight	2 kg with battery
Dimension	175 x 225 mm
Pan	SS 304
Body	ABS Mould
Vision	Audio & Visual Alarm
Volume	Set volume

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country ;
Distributors & Channel
partners World Wide