

CR - 3006

COMPUTED RADIOGRAPHY




EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► TECHNICAL SPECIFICATIONS



Cassette	CR Cassette 14 x 17 (inch) CR Cassette 10 x 12 (inch) CR Cassette 18 x 24 (cm)
Slots	Single Cassette feed
Dimensions (W x D x H)	436 x 636 x 196mm
Weight	21.5kg (47.4lbs)
Grayscale Resolution	Acquisition : 16 bits per pixel Display : 12 bit per pixel
Power Supply Conditions Network PC Connection Computer Min. Requirements	Single Phase 50 ~ 60 Hz AC 90 ~ 264V AC 90 ~ 264V 100 M Bit USB 2.0 OS Win 7 or 10 CPU : Intel i3 RAM : 4 GB HDD : 500 GB Monitor : More than 17" LED Monitor
Operating Conditions	Temperature : 10 ~ 30°C Humidity : 15 ~ 82% RH

►► Overview

We offers such a better clinic environment with no chemicals, ideal space, high resolution image quality and affordability.

Compact & Lightweight

- The length is just 636mm, so you can put it on every normal desk
- It weight only 21.5kg (47.4 lbs)

Durability

- Extreemely simple structure design
- Cassette sorter protects the inside of scanner from dust

Authentic Rigid Type Among Clinic

- No damage or scratches to the image plate during scanning & erasing
- Scan & Erase without roller
- No cut-off image in winter or cloud conditions

Cassette

Strong Structure

- Strong enough against external impact
- Able to bear up to 150kg of weight on cassette

Featherlight


- 14 x 17" : 2.1kg / 10 x 12" : 1.3kg / 18 x 24cm : 0.067kg

Dust free & Easy Cleaning

- Easy to clean up dust on Ips
- Prevent dust from outside

Barcode System

- Automatic cassette size recognition and scanning



User Friendly Design

- Barcode label
- Various colors
14x17(Green), 10x12(Blue)
18x24cm(Pink)



►► Overview (Software)

About our software, is simple, fast, stable and comootable.

Emergency study

- Instantly scan in worklist without creating patient information
- Drag & Drop scanned image to patient list



Easy multilingual support

- It can be easily applied by changing the english list to the corresponding language



Intuitive UI



Lightweight and Robust

- Runs naturally on low - end computers and monitors

Preview

- Can see real-time scanned image
- Decide immediately to repeat

Compatibility

- Can connect with any DICOM printer and paper printer
- Up to 4 division possible



Simplicity

- Easily understand the shooting direction by displaying the tube label as the shape of the cassette plate

Plenty of editable Tools



- Flip / Rotate
- Distance
- Angle
- Rectangle
- Marketing
- Pan / Zoom
- Simple image processing

Special Tools

- Rotate to the desire angle
- Auto Stitching (Coming Soon)
- Client Viewer (Coming Soon)



HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies




UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Partical
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide