

CC - 3312

Cyto Centrifuge



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► PRODUCT INTRODUCTION

It is controlled by microprocessor. equipped with DC Brushless motor. digital and low noise. It can be used to take erythrocyte-serum test, identify antigens and antibodies, and give Other relative tests in immune hematology examination and research laboratories.

►► FEATURES

- Microcomputer control, brushless OC motor. running low noise. high speed accuracy.
- Automatic balance, with large capacity centrifugal cavity, the temperature rise is small,
- Digital display, user-friendly interface, the operation more convenient,
- Equipped With electronic door locks, door cover protection, speeding, unbalance and other protection; fault automatic alarm recognition. safe and reliable,
- 5 minutes can be 8-12 specimens, simple operation, no pollution, cell smear uniform, clear background,

►► TECHNICAL DATA

Model No	CC 3308	CC 3312
Max: Speed	2200rpm	2200rpm
Max. RCF	890 xg	890 xg
Max, Capacity	8 Sample	12 sample
Cytofunnel Height	50 mm	50 mm
Timer	1-99 min	1-99 nm
Noise	<55db	<55db
Motor	Converter motor	Converter motor
Motor Power	200W	200W
voltage	AC 220V, 50/60Hz, 5A	AC 220V, 50/60Hz, 5A
N.W.	22 kg	22 kg
G.W	26 kg	26 kg
Machine Dimension	540X440X300mm (LxWxH)	540X440X300mm (LxWxH)
Package Dimension	620x520x390mm (LxWxH)	620x520x390mm (LxWxH)

►► ROTOR

 <p>No. 1</p>	<p>Max speed : 2200r/min Max capacity : 12 PCS CuPS Max rcf : 890xg</p>	
 <p>No. 2</p>	<p>Max speed : 2200r/min Only for CC 3308</p>	
 <p>No. 3</p>	<p>Max speed : 2200r/min Only for CC 3312</p>	

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies




UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Particle
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide