

# Deep Freezer

## DP-2388

### SPECIFICATION:

- Capacity: 400 Litres
- Dimensions: 1270W 620D 845H mm
- Input Power: 1800 W
- Weight: 55kg
- Temperature range: -10°C to -86°C

### SPECIAL FEATURES:

- Upright type.
- Inside temperature ranging from -10°C - 86°C can be set freely controlling precision is 1oC.
- Stainless steel interior, Outer is painted steel panel
- Front opening lockable door with full height handle.
- 155 mm foaming insulation.
- Vacuum release port
- Inner 3 shelves stainless steel with inner door.
- Imported Danfoss compressor, Two compressors parallel connection to reduce noise and improve the cooling temperature come to -86°C steady
- 2-4 hour come to -86°C very fast
- Fan electromotor with internationally famous Germany EBM.
- Microprocessor-based temperature control system, platinum resistance sensors.
- Digital temperature display.
- Keyboard lock and password protected configuration page.
- Audible and visual alarm:
- High or low temperature alarm, Power failure alarm, Low battery alarm Door ajar alarm, Filter blocking alarm, system failure alarm.
- With temperature printer, record 7 days data.
- 72 hour battery back up for printer and controller display.
- 4 units Casters for easy handling.
- Mixture gas with R134a, R404a, R23, R24, N2, R600a and more CFC and HCFC free.
- Certificate: ISO9001, ISO14001, ISO13485.
- Power supply: 220V/50 Hz 1 phase, can change as 220V 60Hz or 110V 50Hz or 110V 60Hz
- Optional: Chart Recorder, CO2 Backup, Storage racks, boxes.



## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



# Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [info@analyticalfoundation.org](mailto:info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## Reach us @



HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W. [www.ais-india.com](http://www.ais-india.com)  
[www.analycalgroup.net](http://www.analycalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country ;  
Distributors & Channel  
partners World Wide