


# DDA - 3200

## Digital Diagnostic Audiometer



## ►► Product Introduction

Digital Diagnostic Audiometer is advanced micro-controlled audiometer incorporating the latest in digital technology. Digital Diagnostic Audiometer is a compact portable designed with Advanced Digital Technology. It generated precise, distortion free tones and accurate attenuation steps for speedy, genuine & accurate testing procedures.

## ►► Product Introduction

Features	
Test	Pure Tone (Air Conduction, Bone Conduction)
Masking	Wide Band Masking
Speech Test	Speech Test (-10dB to 120dB)
Special Test	Pulse Tone, Tone Decay
Talk Over	Talk Over Amplifier Inbuilt
Calibration	Individual Calibration of Phone L, Phone R & Bone
Distortion	Tone : <0.1%
Key Functions	Feather Touch Control soft keys
Display	Backlite LCD Display, Two Rows of 16 Characters (16 x 2)
USB Port	PC Connectivity through USB Port
Dimensions	200 (L) X 180 (W) X 70 (H) mm
Weight	1 Kg.
Power Supply	AC/DC Adaptor.200-230V 12V-2Amp.DC (Output)

## ►► Digital Diagnostic Audiometer

Air Conduction	Bone Conduction
125 Hz: - 10dB to 80dB	125 Hz: - 10dB to 50dB
250 Hz: - 10dB to 90dB	250 Hz: - 10dB to 50dB
500 Hz: - 10dB to 120dB	500 Hz: - 10dB to 80dB
750 Hz: - 10dB to 120dB	750 Hz: - 10dB to 80dB
1 KHz: - 10dB to 120dB	1 KHz: - 10dB to 80dB
1.5 KHz: - 10dB to 120dB	1.5 KHz: - 10dB to 80dB
2 KHz: - 10dB to 120dB	2 KHz: - 10dB to 80dB
3 KHz: - 10dB to 120dB	3 KHz: - 10dB to 80dB
4 KHz: - 10dB to 120dB	4 KHz: - 10dB to 80dB
6 KHz: - 10dB to 110dB	6 KHz: - 10dB to 80dB
8 KHz: - 10dB to 100dB	8 KHz: - 10dB to 80dB

## HPLC Servicing, Validation, Trainings and Preventive Maintenance :

**HPLC Servicing** :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

**Trainings** :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

**AMC's/CMC** :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

**Validations** :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

**Instruments** :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



## About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

## Our Products & Technologies



UV/VIS  
Spectro 2080+  
Double Beam



Infra FTIR



Optima Gas  
Chromatograph  
3007



Optima Gas  
Chromatograph  
2979 Plus



Flash  
Chromatograph



Atomic Absorption  
Spectrophotometer



Liquid Particle  
Counter



Optical Emission  
Spectrophotometer



DSC/TGA



Semi Auto Bio  
Chemistry Analyzer



HEMA 2062  
Hematology  
Analyzer



Micro Plate  
Reader/Washer



URINOVA 2800  
Urine Analyzer



Total Organic  
Carbon 3800



Fully Automated  
CLIA



NOVA-2100  
Chemistry Analyzer



PCR/Gradient PCR/  
RTPCR



TOC  
Analyzer



Laser Particle  
Size Analyzer



Ion Chromatograph



Water purification  
system

## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



**Analytical**  
**Foundation**

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [Info@analyticalfoundation.org](mailto:Info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.



 **Analytical**®  
Technologies Limited

**HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)**

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)

W. [www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide