

# Double Weighing Balance DWB-3006



## Features:

Operatore put two buckets (with blood bags in each bucket) on both plates

Display shows weight difference and an indication lamp at lighter side turns on.

Operators just add a few weights - chips (e. g. 1 gram each) into bucket until GOOD lamp turns on. The lamp turns on when weight difference becomes less than allowance

Accurate simple and very easy to use

Password Protection

## Specification:

- Range of weight measurement: 0-3000 grams/ml
- Weight difference to be displayed: 0-3000 grams/ml
- Weight Difference for GOOD lamp" 2 gram (may vary upon request)
- Weight display panel: GOOD, A, B
- Weight sensor: land cell
- Control: Power ZERO
- Accuracy:  $\pm 2$  gram
- Power consumption: 5 W
- Weighing Range: 2500 grams/ml
- Storage: -20 - 70°C (-4 - 150°F)
- Temp. operation: 0 - 40°C (32 - 140°F)

## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



**Analytical**  
**Foundation**

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [info@analyticalfoundation.org](mailto:info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## Reach us @



<b>Corporate &amp; Regd. Office:</b> Analytical House, # E67 & E68, Ravi Park, Vasna Road, Baroda, Gujarat 390 015. INDIA	<b>HPLC Solutions</b> T: +91 265 2253620 +91 265 2252839 +91 265 2252370 F: +91 265 2254395	<b>MultipleLabs</b>	<b>Analytical Bio-Med</b> E: <a href="mailto:info@hplctechnologies.com">info@hplctechnologies.com</a> <a href="mailto:info@multiplelabs.com">info@multiplelabs.com</a> <a href="mailto:info@analyticalgroup.net">info@analyticalgroup.net</a> <a href="mailto:info@analyticalbiomed.com">info@analyticalbiomed.com</a>	<b>Analytical Distributors</b> W: <a href="http://www.ais-india.com">www.ais-india.com</a> <a href="http://www.analycalgroup.net">www.analycalgroup.net</a> <a href="http://www.hplctechnologies.com">www.hplctechnologies.com</a> <a href="http://www.multiplelabs.com">www.multiplelabs.com</a>	<b>Analytical Foundation (Trust)</b> Sales & Support Offices: across the country : Distributors & Channel partners World Wide
--	---	---------------------	--	---	---