

# M-DYC-3000 Series

## Dry Cabinet



M-DYC-3060



M-DYC-3080



M-DYC-3120



M-DYC-3140



M-DYC-3180



M-DYC-3160



M-DYC-3240-3540



M-DYC-3320-3435

EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

**Analytical Technologies Limited**

An ISO 9001 Certified Company

[www.analyticalgroup.net](http://www.analyticalgroup.net)

ATL Dry Cabinets with thermoelectric Dehumidifying systems are used for storing photographic and optic lenses, film, camera equipment, audiovisual equipment, paper materials, musical instruments, chemical and medical materials, PCB and IC electronic components.

►► **Features :**

- Cold rolled steel plate construction, with 3.2mm Tempered glass viewing windows and key lockable magnetic closing doors.
- LED display with microprocessor control.
- Desiccant dehumidifying system, alternately absorbs and emits moisture, controlled by ICtimer.

Model	Humidity Range (%RH)	Capacity(L)	Shelf(pcs)	Inner Size WxDxH (mm)	Outer Size WxDxH (mm)
M-DYCM-3060	20~60	65	1	398x375x400	400x400x475
M-DYCD-3060	10~20				
M-DYCDE-3060	1~10				
M-DYCM-3080	20~60	84	2	398x375x545	400x400x620
M-DYCD-3080	10~20				
M-DYCDE-3080	1~10				
M-DYCM-3100	20~60	115	2	398x375x685	400x400x760
M-DYCD-3100	10~20				
M-DYCDE-3100	1~10				
M-DYCM-3120	20~60	128	3	398x375x790	400x400x865
M-DYCD-3120	10~20				
M-DYCDE-3120	1~10				
M-DYCM-3140	20~60	147	4	398x375x940	400x400x1015
M-DYCD-3140	10~20				
M-DYCDE-3140	1~10				
M-DYCM-3180	20~60	189	5	398x375x940	400x400x1315
M-DYCD-3180	10~20				
M-DYCDE-3160	1~10				
M-DYCM-3160	20~60	160	3	448x422x848	450x450x1030
M-DYCD-3160	10~20				
M-DYCDE-3160	1~10				
M-DYCM-3240	20~60	240	3	598x370x1083	600x400x1265
M-DYCD-3240	10~20				
M-DYCDE-3240	1~10				
M-DYCM-3320	20~60	320	3	898x420x848	900x450x1030
M-DYCD-3320	10~20				
M-DYCDE-3320	1~10				
M-DYCM-3435	20~60	435	3	898x570x848	900x600x1030
M-DYCD-3435	10~20				
M-DYCDE-3435	1~10				
M-DYCM-3540	20~60	540	3	598x720x1318	600x750x1500
M-DYCD-3540	10~20				
M-DYCDE-3540	1~10				

## HPLC Servicing, Validation, Trainings and Preventive Maintenance :

**HPLC Servicing** :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

**Trainings** :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

**AMC's/CMC** :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

**Validations** :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

**Instruments** :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



## About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

## Our Products & Technologies



UV/VIS  
Spectro 2080+  
Double Beam



Infra FTIR



Optima Gas  
Chromatograph  
3007



Optima Gas  
Chromatograph  
2979 Plus



Flash  
Chromatograph



Atomic Absorption  
Spectrophotometer



Liquid Particle  
Counter



Optical Emission  
Spectrophotometer



DSC/TGA



Semi Auto Bio  
Chemistry Analyzer



HEMA 2062  
Hematology  
Analyzer



Micro Plate  
Reader/Washer



URINOVA 2800  
Urine Analyzer



Total Organic  
Carbon 3800



Fully Automated  
CLIA



NOVA-2100  
Chemistry Analyzer



PCR/Gradient PCR/  
RTPCR



TOC  
Analyzer



Laser Particle  
Size Analyzer



Ion Chromatograph



Water purification  
system

## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



# Analytical Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [Info@analyticalfoundation.org](mailto:Info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## Reach us @



 **Analytical**®  
Technologies Limited

HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W. [www.analycalgroup.net](http://www.analycalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)  
[www.ais-india.com](http://www.ais-india.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide