

DI - 3460

DRY IMAGER



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► TECHNICAL SPECIFICATIONS

The Dry Imager is a thermo-graphic film processor designed to copy and send images through DICOM network protocol. It uses the latest direct dry thermal imaging technology which accommodates for a full range of modalities including CT, MRI, DR, CR, Digital Gastrointestinal, Nuclear Medicine, Mobile X-Ray Imaging, and Dentistry, etc. The Dry Imager dedicates to precision in diagnosis with its outstanding image quality, and offers affordable imaging catering to your needs.

Model	DI-3430	DI-3460	DI-3760
Print Technology	Direct thermal (dry, daylight-load film)		
Spatial Resolution	320dpi (12.6 pixels/mm)		508dpi (20 pixels/mm)
Throughput	14"X17" w"70 sheets/h; 8"X10" w"90 sheets/h		14"X17"≥60 sheets/h; 8"X10"≥80 sheets/h;
Grayscale Resolution	14 bits		
Film Transferring Method	Rubbing		
Film Tray	One supply tray, 100-sheet capacity total	Two supply trays, 200-sheet capacity total	
Film Sizes	8"X10", 10"X12", 11"X14", 14"X17"		
Applicable Film	Medical Dry Thermal Film (blue or clear base)		
Interface	10/100/1000 Base-T Ethernet (RJ-45)		
Network Protocols	Standard DICOM 3.0 connection		
Image Quality	Automatic calibration using built-in densitometer		
Control Panel	Touch Screen, Online Display, Alert, Fault and Active		
Power Supply	100-240V AC 50/60Hz 600W		
Weight	50Kg		
Operating Temperature	5°C-35°C		
Operating Humidity	35%-85%		
Storage Humidity	30%-95%		
Storage Temperature	-22°C-50°C		
Base Holding	Optional		

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Partical
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.



 **Analytical**[®]
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide