

ESR-2020




Easy & Conomical

Plug in & start testing



Smart & Safe

Load specimem



Rapid & Reliable

Results in 30 minutes

- Westergren manual ESR results vary with temperature and humidity.
- Varied Westergren manual ESR results may be obtained with blood sample 'X' if performed at 17°C; 25°C; 32°C and soon.
- With Automated ESR as the westergren ESR results are auto- corrected to 18°C (Manley) in 30 minutes, the above referred sample 'X' will give same result even if the test is done at 17°C; 25°C; 32°C and so on (if performance in accordance with the intructions given by the manufacturer).

Avoid giving wrong Westergren manual ESR results, switch to Automated ESR 2010 and give auto-corrected westergren results at a standard 18°C, irrespective of your laboratory temperature.

Not just Automation, but also Auto-correction of ESR Results

ESR 2020 Automated erythrocyte Sedimntation rate Analyzer

Features

- Automated ESR analyzer based on the modified Westergren method
- ESR results suto-corrected at 18°C
- Reliable and accurate results in 30 minutes
- Excellent correlation to the standard 1 hour Westegren method
- Useful for creating safe and Healthy Environment for the user
- Maintance free

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide