


Electrophoresis Unit 3000Series



►►► Specification:

<ul style="list-style-type: none"> > Dimension: > Number of sample: > Buffer Volume: > Number of casting Tray: > Number of Combs: > Comb Capacity (5mm Deep Well): > Connecting Cord: > No. of Platinum electrodes: > Buffer Chamber with Safety Lid: 	<p>25.5cm(L)x15.5(B)x7cm(H) Upto 36 450ml 7x5cm, 13x6cm, 13x10cm 8-wells 1.5mm thickness comb (Well Width-5mm, Interspace-3mm) 13-wells 1.5mm thickness comb (Well Width-5mm, Interspace-5mm) 18-wells 1.5mm thickness comb (Well Width-4mm, Interspace-3mm) 37l, 132l, 30 l Red and Black (Each One) Positive and Negative (Each One) 1 No.</p>
--	--

►►► Power Supply Specification:

<ul style="list-style-type: none"> > Output Range: > Type of output: > Output terminals\: > Timer: > Start /Stop Function: > Display: > Operating conditions: > Safety features: > Input protection: > Input power (actual): 	<p>10–500 V, fully adjustable in 1 V steps, 1-500 mA fully adjustable in 1 mA steps Constant voltage or constant current 4 pair recessed banana jacks in parallel 1 min–999 min , fully adjustable Yes 3-digit LED 0–40°C, 0–90% humidity in the absence of condensation No-load detection, sudden load change detection, overload/short-circuit detection Fuse on both hot and neutral 220 VAC/50hz</p>
--	---



▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

▶▶▶ Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.