

EW - 3000

Embedding Workstation
with biospecimen traceability



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

►► EW- 3000 - EMBEDDING WORKSTATION - MAIN FEATURES



- Holds 5 liters of molten paraffin
- Programmable heating
- Reservoir & conduits always heated, ready for urgent operation
- Dispenser compatible with Mega and SuperMega cassettes



- Holds up to 300 cassettes
- Milestone racks and any standard rack compatibility
- Two built-in paraffin trimmers



- Membrane keys
- Programmable daily heating hours including non working holidays
- Extra touch screen tablet option with bar code reader for case traceability



Storage Capacity of up to 400 molds
Programmable heating
Space suitable for Standard, Mega or
SuperMega cassettes

►► **EW- 3000 - CM - COLD EMBEDDING MODULE - FEATURES**

- Large cooling surface
- Programmable temperature
- Consistent temperature distribution over the entire surface
- Stand-alone operation possible
- Can be placed either on the right or the left side of the dispensing console



Temperature from 0° down to -12° Celsius
Space for 60 standard or 45 Synergy molds
Environmentally safe refrigerant (R134a)



►► EW- 3000 - EMBEDDING WORKSTATION - MAIN INCLUDES

- EW- 3000 Dispensing Module
- EW- 3000 Thermal Module
- EW- 3000 Cold Embedding Module*

►► ENHANCE SPECIMEN TRACEABILITY VISUALLY



Options

- HistoDream PC: Tracking system terminal, including Viewer Software and 2D bar code reader
- Magnifying glass
- Heated tweezers 1 mm, 2mm, 4mm

Technical Specifications

- Total Dimensions: 102 cm X 64.1 cm X 40.5 cm (WxDXH); (40.16" X 25.24" X 15.95")
- Total Weight: 61 Kg (134.5 lbs)

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies




UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Particle
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide