

FUME HOOD MODEL:3000 SERIES

FH series Fume hoods are the first defense to minimize chemical exposure to research workers. They are considered the primary means of protection from inhalation of hazardous vapors. Our FH series have three models in total, the details, please see the following,

▶▶▶ Air Pattern and Protected Area:

Advantage:

1. UV lamp for sterilization.
2. LCD displayscreen. See the air velocity directly.
3. The fan is inside. Easy to install, save time.
4. The air velocity is adjustable (6 levels).
5. Electrical controlled glass.
6. Removable work surface. Easy to clean.
7. Standard accessories: Taps, sink, duct, sockets, UV lamp, fluorescent lamp, base cabinet.
8. Working bench : 18mm Jet black granite platform with under bench storage cabinet



▶▶▶ Specification:

Model	3000 Series	
External Size (W*D*H)	1000*750*2200mm	
Internal Size (W*D*H)	820*620*750mm	
Work Surface Height	620mm the exterior size include the cabinet height	
Protection Class	Class I	
Protection Type	IP 20	
Air Velocity	≥0.6m/s	
Noise	≤60dB	
Exhusted Duct	PVC standard length, 4 meter	
Pipe collar	Φ290	
Pipe Strap	1 Pc	
Fluorescent Lamp	10W*1	
UV Lamp	20W*1	
Blower	centrifugal blower, Speed Adjustable, the fan inside fume hood	
Glass Window	Two layer toughened glass, 5mm thickness, motor control height adjustable.	
Power Supply	110-220V/50-60	
Consumption	800W/1400W	
Material	Exterior body: Cold roll steel Sheets, coated with antibacterial powder, interior, Adopts high grade melamine board with good acid and alkali, resistance function work table using solid physiochemical board	
Standard Accessory	Water tap, Air tap, Water Cup, Base cabinet	
	Waterproof socket*2	
Optional Accessory	HEPA filter	
Net Weight	200kg	
Gross Weight	260kg	
Package Size (W*D*H)	Main Body	1.15*0.9*1.72
	Base Cabinet	1.14*0.9*0.9

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide