

HPLC-P6000 Series

HIGH PRECISION DOUBLE PISTON PUMP



►► INTRODUCTION

P6000 HPLC pump is a precision delivery system that working chiefly by means of CPLD technology, and controlled by micro-processor. It based on all-digital frequency generating system, advanced motor driver system and friendly man-machine interface. Operating and controlling become more and more easily. Advanced pump head and check valve combine the most successes experience in the world today and private user's high performance of delivery precision and reproducibility.

EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► P6000 HPLC Pump working mode

P6000 HPLC pump can be configured depend on what the users want:

Iso-Pump for preparative use	P6100 0.1~99.9ml/min, P6250 0.1~249.9ml/min P6500 0.1~499.9ml/min, P61000 0.1~999.9ml/min
Used as HPG Mode for preparative application	In combination with additional P6000 pumps, system capability can be expanded up to quaternary high-pressure gradient.

►► FEATURES

Interchangeable Pump head	Provide 100ml, 250ml, 500ml and 1000ml pump head according to the using environment
Washing Piston rod	The pump working pressure can be automatic calibrated only by press a key before using.
Interactive interface	Using 320x240 LCD display, it is very easy for user to edit programs and more information will be appeared in the screen
Product's information	By checking the product's information, you can find the serial Number, the CMC permission number, the installation date, the total running time and etc. It helps the user decide when the instrument should be maintenance.

►► SPECIFICATIONS

Delivery method	double pistons,main&auxiliary
Maximum working pressure	P6100 20Mpa, P6250 15Mpa, P6500 10Mpa, P61000 5Mpa
Flow Range	P6100 0.1-99.9ml/min, P6250 0.1-249.9ml/min P6500 0.14-99.9ml/min, P61000 0.1-999.9ml/min
Flow rate accuracy	±2%
Flow rate precision	RSD<0.1%
Display	320*240 blue background color
Protection	Pmix&Pmin Data automatically saved after power down
System connection	1/8 capillaries
Control	Stand-alone mode by front keyboard HPG mode:one pump works as master,other works as slave RS232 by PC using communication protocol
Dimensions (W x H x D)	400x270x1800mm
Weight	15kg

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.


Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Particle
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.