

Horizontal Rotary Shaker



> Specification:

Platform Size	24" x 24"
Max. Volume	16 ltrs.
Power Supply	240V, 50Hz
Body	Mild steel platform, painted & with neoprene rubber discs for holding flasks.
Motor	Variable speed motor
Speed Control	Digital rpm controller (20 to 300 RPM)
Capacity	150ml. 49Nos
Timer	Automatic digital timer for control of shaking time.
Features	It accommodates 50ml flask/100ml flask /250ml or 500ml flask.
Clamps for holding Conical flask	

➤ Regulatory compliances



➤ Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

➤ Reach us @



Corporate & Regd. Office: Analytical House, # E67 & E68, Ravi Park, Vasna Road, Baroda, Gujarat 390 015. INDIA	HPLC Solutions T: +91 265 2253620 +91 265 2252839 +91 265 2252370 F: +91 265 2254395	MultipleLabs	Analytical Bio-Med E: info@hplctechnologies.com info@multiplelabs.com info@analyticalgroup.net info@analyticalbiomed.com	Analytical Distributors W. www.ais-india.com www.analycalgroup.net www.hplctechnologies.com www.multiplelabs.com	Analytical Foundation (Trust) Sales & Support Offices: across the country : Distributors & Channel partners World Wide
--	---	---------------------	--	---	---