

IEB - E3154

ICU Electric Bad



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Function Description :




Back-Section



Foot-Section



Trendelenburg



Reverse-Trendelenburg



HI-LOW

IEB-E3154 is an economical five-function electric hospital ICU bed , The remote control controls all movements, It is a recovery platform that can perfectly serve as the hospital' s ICU bed with the below excellent features.

- Four split side rail, drop down under platform.
- Hand remote control.
- Detachable head/foot board with easy lock.
- Caster diameter of 125mm with brakes.
- Bumper caster on four corners.
- Auto-counter to provide the patient from sliding onto foot end.
- With angle indicator embedded in the railing for back and Trendelenburg adjustment (optional).
- Detachable head foot board for emergency operation.
- Manual CPR on both sides(optional).

01 Side rail lock

Tuck-away side railing meets safety standard, special protection mechanism can prevent patient to drop down the railing from inside bed.



02 Remote control

The remote control with clear icon is beneficial for medical staff which can control the bed position from a remote distance.



03 Angle Indicator

Angle indicator on the guardrail. backrest railing indicator for back section adjustment with visible 30 degree mark, knee railing indicator to show degree of trendelenburg.



04 IV hole

There are infusion holes at the head and foot of the bed, so the nurse can infuse in two positions at the same time, which is more convenient.



05 Castor brake

All four castors are equipped with brakes, which are stable and durable.



06 Anti-collision angle

There are anti-collision angles on the four corners of the bed to reduce the impact of hitting the wall.



►► Technical data

Size	2240*1060 *500-780mm
Lying surface	2045* 900mm
High-Low	500-780mm
Back section adjustable	0° - 75°
Knee section adjustable	0° - 40°
Trendelenburg	0° - 12°
Reverse trendelenburg	0° - 12°
Safe working load	170kg
Max working load	225kg

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.


Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Partical
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**[®]
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.