

BED - 3151D

ICU Hospital Bed



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Automatic Hospital ICU Bed



BED-3151D electric hospital bed is full electric operation type ICU bed to take care of the patient in ICU area. It will provide convenient support for both caregiver and patient. The road to recovery is often a critical and sensitive experience. Both patient and medical staff will have to work hand-in-hand to achieve their goal of nursing the patient back to health.

- Motor Quiet and robust electric actuators provide reliable operation
- Build-in controller can control all electrically operated functions and movements
- Electric CPR function control by one single button on the side rails 4 ABS side rails with Tuck-away model, collapsible below the bed surface
- Nurse control located in the footboard. with side rail control and a functions lockout feature
- Auto-counter to provide the patient from sliding onto foot end
- Bed Frame Electrophoresis and powder coating
- Bed platform removable ABS platform
- Head and Footboard ABS material, detachable Central brake system with a wheel diameter of 125mm
- Hand remote control
- Patient control inner constructed in side rails
- Manual CPR

NO.	ITEM	PARAMETERS
01	Overall size	2200*980 *450-750mm
02	Lying size	1950*850 mm
03	Back section up/down	0-70°
04	Foot section up/down	0-40°
05	Height adjustment	450-750mm
06	Trendelenburg	0-12°
07	Reverse-Trendelenburg	0-12°
08	Loading capacity	220kg





01 Bed board : The bed surface is detachable and has a handle for easy cleaning also it with vent holes for ventilation.

02 Side rail lock : Tuck-away full length railing meets safety standard, special protection mechanism can prevent patient to drop down the railing from inside bed.

03 Siderail control : All functions can be operated by the built-in side rail control buttons, which are located inside and outside of the side rail, which is convenient for patients and medical staff to operate.

04 Foot board control panel : Nurse control panel embedded onto the foot board, railing control operation can be locked accordingly when position operation from patient is not wanted.

05 Detachable head and foot board : The detachable head and foot of the bed are equipped with safety locks to lock the head and foot board when it is not in use, and it can be quickly disassembled, which is convenient for medical staff to operate.

06 Central lock : Central lock system is convenient for the nurse to push the patient to the required position and quickly lock it with one key, which can appropriately reduce the workload of the medical staff.

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Particle
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.