

# IR Thermo-3306

Non Contact Forehead IR Thermometer



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

**Analytical Technologies Limited**

An ISO 9001 Certified Company

[www.analyticalgroup.net](http://www.analyticalgroup.net)

## ►► FEATURES:

The forehead gun (infrared thermometer) is designed for measuring the forehead temperature of the human body. It is very simple and convenient to use. Accurate temperature measurement in 1 second, no laser spot, avoiding potential damage to eyes, no need to contact human skin, avoid cross infection, one-click temperature measurement, and check for flu. Suitable for home users, hotels, libraries, large enterprises and institutions, can also be used in hospitals, schools, customs, airports and other comprehensive places, and can also be provided to medical staff in the clinic.

- English system
- Support Fahrenheit and Celsius switching
- Measure Body and Surface temperature
- Memorization of the last 34 measurements
- Three color backlighting (below 37.5°C, show green, 37.5~38°C, show orange; above 38°C, show red)
- Non-contact measuring temperature accurately.
- It can reduce the risk of virus spread effectively.
- Reading temperature display quickly. The response time is only 0.5s.
- Accuracy is more or less than 0.3Deg.C (6Deg.F).
- 34 Measurement data can be saved.
- Different temperature range and different back light.
- Safe to use: When the temperature is over than a certain degree, this thermometer will ring at once.

## ►► SPECIFICATIONS:

- Accuracy:  $\pm 0.3^{\circ}\text{C}$  ( $0.6^{\circ}\text{F}$ )
- Dimension: 145×80×40mm
- Body temperature range: 32.0 ~ 43.0°C (90~ 109°F)
- Operating temperature : 16°C~40°C(61°F~104°F)
- Surface temperature accuracy:  $\pm 0.1^{\circ}\text{C}$  ( $\pm 0.1^{\circ}\text{F}$ )
- Surface temperature range: 0~100°C(32°F~212°F)
- Resolution: 0.1°C
- Response time: 0.5s
- Measuring distance: 5~8cm
- Auto power off : 7s
- Battery: 1.5v of 2 sections

## HPLC Servicing, Validation, Trainings and Preventive Maintenance :

**HPLC Servicing** :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

**Trainings** :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

**AMC's/CMC** :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

**Validations** :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

**Instruments** :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



## About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

## Our Products & Technologies



UV/VIS  
Spectro 2080+  
Double Beam



Infra FTIR



Optima Gas  
Chromatograph  
3007



Optima Gas  
Chromatograph  
2979 Plus



Flash  
Chromatograph



Atomic Absorption  
Spectrophotometer



Liquid Particle  
Counter



Optical Emission  
Spectrophotometer



DSC/TGA



Semi Auto Bio  
Chemistry Analyzer



HEMA 2062  
Hematology  
Analyzer



Micro Plate  
Reader/Washer



URINOVA 2800  
Urine Analyzer



Total Organic  
Carbon 3800



Fully Automated  
CLIA



NOVA-2100  
Chemistry Analyzer



PCR/Gradient PCR/  
RTPCR



TOC  
Analyzer



Laser Particle  
Size Analyzer



Ion Chromatograph



Water purification  
system

## Regulatory compliances



## Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



**Analytical**  
**Foundation**

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [Info@analyticalfoundation.org](mailto:Info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## Reach us @

