

IR Thermo-3306i

Infrared Forehead Thermometer



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

- Adopts reliable sensor
- New & patented probe design ensures high accuracy
- Excellent adaption to ambient temperature
- Accurate and reliable even under complicated surrounding
- Body mode and object mode available. it can also can measurement amient temperature.
- Fever alarm
- Large LCD with 3 colors backlight
- Automatic shut off.

Range

- Human body : 32.0°C - 42.9°C(89.6°F-109.2°F)
- Object : 0°C-100°C(32°F-212°F)

Accuracy

- Humana body : $\pm 0.2^{\circ}\text{C}/0.4^{\circ}\text{F}$
- Object : $\pm 1.0^{\circ}\text{C}/1.8^{\circ}\text{F}$

- **Resolution** : $0.1^{\circ}\text{C}/0.1^{\circ}\text{F}$
- **Working condition** : $16.0^{\circ}\text{C}\sim 35.0^{\circ}\text{C}$ ($60.8^{\circ}\text{F}\sim 95.0^{\circ}\text{F}$) RH \leq 85%
- **Storage condition** : $-20.0^{\circ}\text{C}\sim 55.0^{\circ}\text{C}$ ($-4.0^{\circ}\text{F}\sim 131.0^{\circ}\text{F}$) RH \leq 93%
- **Power supply** : d.c. 3V 2AAAalkaline Battery

Power consumption

- When off \leq 10uW
- When measurement \leq 30mW

- **Power level indicator** : Indication for low power level
- **Memory** : 1 set reading recall
- **Display** : LCD back-lit display
- **Reading scale** : Celsius of Fahrenheit
- **Automatic shut off** : In 20 seconds
- **Dimensions** : 141 mm x 86 mm 39 mm
- **Net weight** : 82 g
- **Service life** : 2 years
- **Measurement method** : Non-contact
- **Effective distance** : 3 cm - 5 cm



Precision



Backlit Display



Quick Measurement



Fever Alarm



Forehead Type



Memory Recall

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @

