

Infrared Thermometer



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

- Dual temperature mode
- 10 set memories
- Three color backlight
- Automatic shut-down
- Accuracy $\pm 0.2^{\circ}\text{C}$



Measuring Accuracy $\pm 0.2^{\circ}\text{C}$



High Fever Alarm



Automatic shut-down

►► Packing Specifications

- Net Weight : 186g
- Host Size : 17.5cm * 12.5cm * 5.2cm(100pcs/corton)
- Gross Weight : 24kg
- Packing Size : 43cm * 52cm * 62cm

►► Specifications

Measuring Unit:

- $^{\circ}\text{C}/^{\circ}\text{F}$

Resolution:

- 0.1°C

Measuring Range:

- Forehead Temperature Mode 32.0 °C ~ 42.9 °C;
- Surface Temperature Mode 0 °C ~ 99.9 °C

Clinical Accuracy:

- Clinical Deviation ± 0.2 °C;
- Clinical Repeatability ± 0.3 °C;

Maximum Permissible Error:

- 32.0°C~34.9°C, ± 0.3 °C;
- 35.0°C~42.0°C, ± 0.2 °C;
- 42.1°C~42.9°C, ± 0.3 °C.

Power Supply Voltage:

- DC3V(2*AA)

Memory Storage:

- 10 sets of latest measurement data

Automatic Shutdown:

- Within 40 seconds

Operating Enviroment:

- Temperature : 16°C ~ 35°C;
- Humidity : $\leq 85\%$
- Atmospheric Pressure : 70kPa ~ 106kPa

Transportation And Storage Enviroment:

- Temperature : ≤ -20 °C ~ 55°C;
- Humidity : 10% ~ 93%
- Atmospheric Pressure : 60kPa ~ 106kPa

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.