

Laboratory Furniture

Fume Hoods

3000 Series



EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

Laboratory Furniture



Laboratory Fume Hoods



Laboratory Accessories

►► Laboratory Furniture



Laboratory Furniture: Perfect for today's lab versatile enough for tomorrow's changes. You can specify work surfaces, shelves and back splashes in a variety of finished surfaces, including Granite, Chemical-Resistant Laminate or Stainless Steel - depending on your lab requirements. Painted or Stainless Steel modular cabinetry of the highest quality fit and finish with hundreds of drawer and door combinations means you can store lab items exactly where you need them.



Island Bench & Work Benchese: Analytical lab Workbenches are designed and manufactured to offer unmatched quality, value, and durability to provide superior performance for today's lab environments. Work Surfaces and shelves can be specified in materials such as High-Pressure Laminate, Stainless Steel, or Composite Resinwood Surfaces. Base Cabinets and Support Legs, used to form the foundation for every Island Bench / Work Bench, are available in Painted Steel or Stainless Steel. Any Analytical lab Base Cabinet can be grounded on : (1) Pedestal Frame (2) C - Frame (3) H - Frame (4) Structure Frame.



►► Laboratory Fume Hoods



Analytical lab offers lab furniture-solutions for almost complete types of laboratory to clarify, if the Furniture will be needed for already existing premises resp. labs, or if it is about a new build. In the planning phase on new construction, by request Analytical lab can impinge on the layout planning, to straighten the intersections for media, water and power optimally on the later workplace design. The same is true for the interior design (Pull in of subceilings, walls). The challenge on already existing labs is, to adapt the new furniture to the existing, on material, shape and colour, so that it comes into being a harmonic impression (compatibility). Important parameters will be obtained off the layout for the lab furniture planning:



- Usage space (length x width)
- Height ceiling
- Load-bearing walls
- Site staircase/ hallway
- Site and measures of doors and windows
- Parapet level window
- Fixed water connections
- Number rooms

►► Laboratory Accessories

Laboratory Furniture / Laboratory Benches

- General Features.
- Modules are for standing & Sitting Postures.
- Standard construction of cold rolled steel finished with durable epoxy polyester powder coating.
- Assemblies are fully modular with flush face construction.
- All drawers run on smooth & durable telescopic 2 stage pull out channels.
- Drawer front are double skin construction.
- Modules on pedestal frames.



►► Laboratory Fume Hoods / Fume Cup-Boards with Extraction System



Fume cupboard represents a low-voltage appliance. Which is necessary for the correct function of this instrument. Since the work in fume cupboard is mainly associated with hazardous chemical agents regular service are recommended.

Fume cupboard: outer dimensions - width 1200, 1500, 1800, depth 750 / 900, height 2250 mm. The worktop is 900 mm. above the floor. The material is made of 1.5 mm G.I. thick steel sheet; surface is epoxy painted according to the RAL scale coating made with resin layer – Front window in the frame is vertically sliding, Toughened glassed vertical or horizontally sliding in the frame .



BENCH TOP



WALK-IN



DISTILLATION

►► Regulatory compliances



►► Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

►► Reach us @

