



**Analytical**<sup>®</sup>

Technologies Limited

## HD LAPAROSCOPE 10mm 0 Degree



### Specification

- Wide Angle Full Screen,
- Optimal center-to-edge resolution for enhanced picture quality.
- Angle of view: 0°
- Diameter 10mm
- Length 33 cm
- Fiber optic light transmission incorporated
- Standard ocular window for coupling the camera head
- Scratch resistance tip lens
- Advanced Rod lens system for optimum brightness, contrast & definition



## 1-CCD CAMERA VM-2008

### Specification

- Smart And User-Friendly
- Single Chip Camera System
- Video Outputs: S-Video and Composite Video
- The Camera offers Signal to Noise ratio range of 60-70dB
- Camera Head offers a focusing coupler for even focus control



## **MEDICAL GRADE MONITOR**

### **Specification**

- Colored Monitor 24"
- Flat Panel Monitor
- PAL system compatible.
- With DVI inputs
- Compact & Lightweight design.
- Resolution more than 1100 lines

**DW  
TOR**

## Specification

- High flow of 30 liters or more with LED display
- Microprocessor controlled & Software driven
- Soft approach pressure control for safe recovery of abdominal pressure
- Have alarms with min 0.1 L flow rate
- Internal leakage detection capability
- Having internal venting system for safety
- Unit include hose



## Irrigation Pump VM15-1733



### Specification

#### Pressure

- Pump Capacity max. : 400 mm Hg.
- Irrigation rate max. 1.8 L / Min (+/- 10%)

#### System

- Suction rate max.1.8 L/ Min (+0.4 / 0.1 L /Min)
- 100 V ~ / 110~150 V~ / 220-240 V- 50/60 Hz

#### Dimension / Weight

- Dimension (W/H/D) : 280 X 137 X 408 MM
- Weight : 5.9 kg (Approx.)

### Features

- Aspirates & Irrigates upto 3 l / min
- Large two litre aspiration container.
- Special filter ensure sterility
- No special adjustments necessary, aspiration & irrigation are regulated and selected by the S/I instruments.
- Integrated overflow protection

## Regulatory compliances



## Corporate Social Responsibility



**Analytical**  
**Foundation**

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [Info@analyticalfoundation.org](mailto:Info@analyticalfoundation.org)
2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.
3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## Reach us @



 **Analytical**®  
Technologies Limited

HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W. [www.ais-india.com](http://www.ais-india.com)  
[www.analycalgroup.net](http://www.analycalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide