

## LEEP SYSTEM



## ATL 2010 LE

## Specification

- A. Electrosurgical generator with isolated power output LED display Located in Front for precise power selection, deliver and easy to use.
- B. Provision of choice to CUT, BLEND and COAG. Wave form to accommodate subtle differences in technique and electrode performance.
- C. RF output frequency 350 - 450 KHz power cut 0-100 watt.
- D. Flash faceplate membrane facilitate operation cleaning.
- E. Microprocessor controlled for precision, accuracy, reproducibility and safety.
- F. Pneumatic / normal spring type pedal for maximum safety.
- G. Audible safety features include distinct tones for each operating setting.
- H. Automatic self test mechanism ensures accurate system operation.
- I. High air flow efficiently captures smoke plume with a variable speed control.
- J. Triple stage filtration captures airborne particulate matter, vapor and odor with a 99.999% efficiency level.
- K. Virtually maintenance free.
- L. Replacement filters available.
- M. Standard accessories are.
  - (1) Hand piece adaptor
  - (2) Patient return (single use)
  - (3) Smoke evacuator package
  - (4) Smoke evacuator pre filter
  - (5) Smoke evacuator reducers
  - (6) Smoke evacuator disposable tubing (6 ft.)
  - (7) Ball electrode
  - (8) Electrode 2cm x 0.8cm, 12c"

## Support & Services

### HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing	:We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.
Trainings	:We also take up preventive Maintenance to reduce downtime of HPLC's
AMC's/CMC	:We offer user training both in-House and at customer sites on HPLC principles, operations, troubleshooting.
Validations	:We have protocols for carrying out periodic Validations as per GLP/GMP/U SFDA norms
Instruments	:We offer instruments / Renting Services Modules like pumps, detector etc. on Rent.



### About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO :9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

### About Analytical Technologies



UV/VIS  
Spectro 2060+  
Double Beam



PTIR



Gas Chromatograph  
3000



Gas Chromatograph  
2979 Plus



Flash  
Chromatograph



Atomic Absorption  
Spectrophotometer



Liquid Partical  
Counter



Optical Emission  
Spectrophotometer



DSC/TGA



NOVA 2020 plus  
Automated Bio Chemistry  
Analyzer



HEMA 2020  
Hematology  
Analyzer



Micro Plate  
Reader/Washer



Water purification  
system



Total Organic  
Carbon



Fully Automated  
CLIA



NOVA Basic Semi -Auto  
Chemistry Analyzer



PCR/Gradient PCR/  
RTPCR



Blood Gas  
Analyzer



Random access Analyzer  
for immunoassay  
Proteins & clinical chemistry



Semen Analyzer



## ▶▶▶ Regulatory compliances



## ▶▶▶ Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



**Analytical**  
**Foundation**

1. Research & Innovation Scientist's awards / QC Professional Award: Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [info@analyticalfoundation.org](mailto:info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## ▶▶▶ Reach us @



 **Analytical**®  
Technologies Limited

HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W. [www.ais-india.com](http://www.ais-india.com)  
[www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.