

MP-3210

MULTI PARA MONITOR



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► General Features :

- Light weight, compact and portable
- High resolution with screen size **12.1"** colour TFT display
- Alarm: Audible and visual Alarm with alarm events review
- **Detachable rechargeable lithium battery**
- 480 hours graphic and tabular of all parameter
- 71 alarms events of all Parameters
- 2400 NIBP measurement Data can be stored and recall
- Battery backup **Approx.90 min**
- Apnea alarm 10-40s
- ECG waveforms of 7 leads display on the same screen
- **32 seconds** full — disclosure waveform review
- Arrhythmia analysis and S-T Segment analysis and its detection
- Data and waveforms color can be stored and recall
- Suitable for **adult, paediatric and neonatal patients**
- Battery backup - **2 hours**
- Dimension = **314 X 145 X 264 mm**
- Weight : **3.8 kg**

►► Technical Specifications :

SpO₂

- Accuracy for PR : **+/-2 bpm**

NIBP

Resolution Pressure : **1mm hg**

ECG

Measurement rang: **0.2 mg-0.2 mg**

Pulse rate

Range : **30-250rpm**

Accuracy for PR+/-2 bpm

Resolution : **1 bpm**

- Measurement Range: 0 ~ 100%
- Accuracy for SpO : $\pm 2\%$ (for 70 ~ 100%) and unspecified for 0 ~ 69%

Pulse Rate

- Range: 30 ~ 250 bpm
- Accuracy for PR : ± 3 bpm

NIBP

- Measurement Range: **10 ~ 270 mmHg**
- Modes: Manual / Auto / Continuous with Oscillometric Method

ECG

- Lead Mode: **3 lead or 5 Lead**
- Lead Selection : I, II, III, avR, avL, avF and V
- Heart Rate Range:
 - » Adult : **35 ~ 300bpm**
 - » Paediatric/Neonatal: **15 ~ 300bpm**
- **ECG Waveform : 2 Channels**
 - » Accuracy : **± 1 bpm** or **$\pm 1\%$**
- ST segment detection:

Measurement range: **0.2 mV ~ 2.0mV**

Respiration

- Accuracy : **+2 rpm**
- Method: **RA-LL impedance**
- Adult : **0 ~ 120 breath/min**
- Paediatric/Neonatal : **0 ~ 150 breathe/mm**

Temperature

- Range : **0 ~ 50°C** • Accuracy : **$\pm 0.1^{\circ}\text{C}$**

**Optional: EtCO₂, IBP, Printer Central monitoring software & cable
(LAN networking or wireLess)**

VGA Output

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.


Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph



Atomic Absorption
Spectrophotometer



Liquid Particle
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**[®]
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.