

# 3000 Series

## Mechanical Stirrer



EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

**Analytical Technologies Limited**

An ISO 9001 Certified Company

[www.analyticalgroup.net](http://www.analyticalgroup.net)

## Mechanical Stirrer 3001

### ►► Specifications

Power	100w
Speed Range	starting up -3000rpm
Timer	0-120min or all the time
Voltage	220v 50hz
Max. Stir Capacity	10L
Stir	Stainless Steel Stir
Weight	9kg
Teflon Stir Bar Are Optional	



## Stable Temperature Electric Stirrer 3002

### ►► Specifications

Power	1060w
Speed Range	starting up -3000rpm
Constant temperature	
Temperature range	RT-100°C
Max. Stir Capacity	10L
Stir	Stainless Steel Stir
Voltage	220v 50hz
Weight	11kg
Teflon Stir Bar Are Optional	



## Speed Measuring Electric Stirrer 3003

### ►► Specifications

Digital display speed	
Power	100w
Speed Range	starting up -3000rpm
Timer	0-120min or all the time
Max. Stir Capacity	10L
Stir	Stainless Steel Stir
Voltage	220v 50hz
Weight	10kg
Teflon Stir Bar Are Optional	



## ►► Servicing, Validation, Trainings and Preventive Maintenance :

Servicing	: We have team of service engineers who can attend to any make of instrument promptly @the most affordable cost.
Trainings	:We also take up preventive maintenance to reduce downtime of instrument's Trainings.
AMC's/CMC	: We offer user training both in-House and at customer sites on instrument principles, operations, troubleshooting.
Validations	: We have protocols for carrying out periodic Validations as per GLP/ GMP/USFDA norms.
Instruments	: We offer instruments/Renting Services Modules like pumps, detector etc. on Rent.



## ►► About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO :9001 Certified Company engaged in Designing, Manufaturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

## ►► Our Products & Technologies



UV/VIS  
Spectro 2080+  
Double Beam



Infra FTIR



Optima Gas  
Chromatograph  
3007



Optima Gas  
Chromatograph  
2979 Plus



Flash  
Chromatograph



Atomic Absorption  
Spectrophotometer



Liquid Partical  
Counter



Optical Emission  
Spectrophotometer



DSC/TGA



Semi Auto Bio  
Chemistry Analyzer



HEMA 2062  
Hematology  
Analyzer



Micro Plate  
Reader/Washer



URINOVA 2800  
Urine Analyzer



Total Organic  
Carbon 3800



Fully Automated  
CLIA



NOVA-2100  
Chemistry Analyzer



PCR/Gradient PCR/  
RTPCR



TOC  
Analyzer



Laser Particle  
Size Analyzer



Ion Chromatograph



Water purification  
system

## ►► Regulatory compliances



## ►► Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



# Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [info@analyticalfoundation.org](mailto:info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## ►► Reach us @



# Analytical<sup>®</sup> Technologies Limited

**HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)**

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W: [www.ais-india.com](http://www.ais-india.com)  
[www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide

Note : Company reserves rights to add/delete/modify the contents /technical specifications of the catalogue without prior notice.