

FM-3910

Metallurgical Microscope



EPCC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Overview

Brand Name:	ANALYTICAL
Model Number:	FM-3910
Theory:	Metallurgical Microscope
Drawtube:	Trinocular
Product name:	Professional Metallurgical Microscope with Bright/Dark field
Magnification:	50X-1000X
Application:	institutes and laboratories
Optical System:	Infinite Color Corrected Optical System
Objective:	5X/10X/20X//50X/100X
Eyepiece:	PL10X/25mm
Nosepiece:	Bright & Dark Quintuple Nosepiece
Stage:	4 inch three layers mechanical stage
Condenser:	Swing-out type achromatic condenser (N.A.0.9)
Illumination:	12V/100W halogen lamp house

►► Item	Specifications	FM-3910
Optical System	Infinite Color Corrected Optical System	•
Viewing Head	Siedentopf trinocular head, Erect image, 30° inclined, interpupillary distance: 50mm~76mm; splitting ratio Eyepiece:Trinocular=100:0 or 0:100	•
	Siedentopf trinocular head, Inverted image, inclined at 30°, interpupillary distance: 50mm~76mm; splitting ratio Eyepiece:Trinocular=100:0 or 20:80 or 0:100	°
Eyepiece	High eyepoint wide field plan eyepiece PL10X/25mm, diopter adjustable	•
	High eyepoint wide field plan eyepiece PL10X/25mm, with reticle, diopter adjustable	°
	High eyepoint wide field plan eyepiece PL10X/26.5mm, diopter adjustable	°
	High eyepoint wide field plan eyepiece PL10X/26.5mm, with reticle, diopter	°

Objective	Bright & Dark Field Semi-Apochromatic Metallurgical Objective 5X, NA=0.15, WD=13.5mm	•
	Bright & Dark Field Semi-Apochromatic Metallurgical Objectives 10X, NA=0.30, WD=9mm	•
	Bright & Dark Field Semi-Apochromatic Metallurgical Objectives 20X, NA=0.50, WD=2.5mm	•
	Bright & Dark Field Semi-Apochromatic Metallurgical Objectives 50X, NA=0.80, WD=1.0mm	•
	Bright & Dark Field Semi-Apochromatic Metallurgical Objectives 100X, NA=0.90, WD=1.0mm	•
	Bright Field Semi-Apochromatic Metallurgical Objective 5X, NA=0.15, WD=19.5mm	◦
	Bright Field Semi-Apochromatic Metallurgical Objectives 10X, NA=0.30, WD=10.9mm	◦
	Bright Field Semi-Apochromatic Metallurgical Objectives 20X, NA=0.50, WD=3.2mm	◦
	Bright Field Semi-Apochromatic Metallurgical Objectives 50X, NA=0.80, WD=1.2mm	◦
	Bright Field Semi-Apochromatic Metallurgical Objectives 100X, NA=0.90, WD=1.0mm	◦
Nosepiece (with DIC slot)	Bright & Dark Quintuple Nosepiece	•
	Bright & Dark Sextuple Nosepiece	◦
	Bright Field Sextuple Nosepiece	◦
	Bright Field Septuple Nosepiece	◦

Frame	<p>Reflected & Transmitted body, low-position coaxial coarse and fine adjustment, coarse adjustment distance: 25mm; fine precision: 0.001mm. With coarse adjustment stop and tightness adjustment.</p> <p>Built-in 100-240V wide voltage transformer, double way power output; intensity adjustable by digital set and reset; switch for reflection and transmission; built-in transmitted filters LBD/ND6/ND25).</p>	•
	<p>Reflected body, coaxial coarse and fine adjustment, coarse adjustment distance: 25mm; fine precision: 0.001mm. With coarse adjustment stop and tightness adjustment.</p>	°
Stage	4 inch three layers mechanical stage with glass plate, moving range: 102mm(Y) *105mm(X)	•
Condenser	Swing-out type achromatic condenser (N.A.0.9)	•
Reflected Illuminator	BD reflected illuminator with iris field diaphragm and aperture diaphragm, central adjustable. With filter slot and polarizing slot. With switch for bright and dark field.	•
Lamp House	12V/100W halogen lamp house, center pre-set	•
Other Accessories	Camera adapter: 0.5X focusing C-mount	•

HPLC Servicing, Validation, Trainings and Preventive Maintenance :

HPLC Servicing :HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.

Trainings :We also take up preventive Maintenance to reduce downtime of HPLC's Trainings.

AMC's/CMC :AMC's/CMC :We offer user training both in-House and at customer sites on HPLC principles, operations, trouble-shooting.

Validations :Validations :We have protocols for carrying out periodic Validations as per GLP/GMP/USFDA norms.

Instruments :Instruments :We offer instruments/Renting Services Modules like pumps,detector etc. on Rent.



About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the global. Analytical Technologies Limited is an ISO:9001 Certified Company engaged in Designing, Manufacturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

Our Products & Technologies



UV/VIS
Spectro 2080+
Double Beam



Infra FTIR



Optima Gas
Chromatograph
3007



Optima Gas
Chromatograph
2979 Plus



Flash
Chromatograph




Atomic Absorption
Spectrophotometer



Liquid Particle
Counter



Optical Emission
Spectrophotometer



DSC/TGA



Semi Auto Bio
Chemistry Analyzer



HEMA 2062
Hematology
Analyzer



Micro Plate
Reader/Washer



URINOVA 2800
Urine Analyzer



Total Organic
Carbon 3800



Fully Automated
CLIA



NOVA-2100
Chemistry Analyzer



PCR/Gradient PCR/
RTPCR



TOC
Analyzer



Laser Particle
Size Analyzer



Ion Chromatograph



Water purification
system

Regulatory compliances



Corporate Social Responsibility

Analytical Foundation is a nonprofit organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards/QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at Info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops/Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds,souls and body by means of yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

Reach us @



 **Analytical**®
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com
www.ais-india.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide