

MICROBIO BIOAEROSOL SAMPLING



The MicroBio MB2 Bioaerosol Sampler has a host of features making it ideal for use in research, clean rooms, pharmaceutical validation and throughout hospitals. Fully validated under the UK DTI VAM programme as a Reference Sampler. Capable of using both 90 mm Petri dishes or 55 mm contact plates and a wide variety of media makes the MicroBio a very low cost instrument to use.

The MicroBio MB2 is an invaluable device for validating cleanrooms to ISO 14698 requirements.



Features

Flow rate 100 litres / minute, factory calibrated

d50 = 1.7 µm for 220 hole head

d50 = 1.35 µm for 400 hole head

Sampling volume 25 to 10,000 litres

6 User programmable pre-set volumes

Programmable delay to start sampling up to 3 hours

Tripod mounting facility

Lightweight (750 g) and portable

Battery powered - rechargeable or alkaline

NiMh (4 x AA) cells and charger supplied

Low noise

Sample over 60,000 litres per charge

Uses 90 mm Petri or 55 mm contact plates

Supplied with the MB2 Bioaerosol Sampler

Petri dish and contact plate spring sets

Option of 220 hole x 1 mm stainless steel or 400 x 0.7 mm anodised aluminium sampling head

Printed operating manual

Padded carry bag

Optional mini tripod

Multi-voltage Ansmann charger with UK, Euro, USA, Australia and 12V in-car plugs and pre-charged cells so the MB2 is ready to use out of the box

Certificate of Calibration



▶▶▶ Regulatory compliances



▶▶▶ Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical
Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

▶▶▶ Reach us @



HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W. www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents / technical specifications of the catalogue without prior notice.