

# Milk Analyzer

## AMIA- 8654



LABS / EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

**Analytical Technologies Limited**

An ISO 9001 Certified Company

[www.analyticalgroup.net](http://www.analyticalgroup.net)

## ►► Specifications:

Model	AMIA- 8654		
Measuring time	Measuring Parameters	Standard complete	Additional options
90 Sec	FAT -0.01% -25% SNF -3% -15% Density -1015 - 1040kg/m <sup>3</sup> Proteins -2% -7%	1. Sample visualization on LD display(green) 2. Possibility to measure 1 sample -90 seconds 3 Adapter 220-110 V 4. Spare silicone tubes 5. Syringe 6. Plastic sample holder 7. Calibrations Cow, Sheep, and UHT milk 8. Operation Manual 9. Cardboard transport case 10. CD -Service Pack	Ph option
			Conductivity
60 Sec	Lactose -0.01% -6% Water content -0% -70% Temp. of milk -1°C-10°C Freezing point - -0.400 -0.700°C Solids - 0.4 -1.5%		High Fat 45% (High Density)

## ►► HPLC Servicing, Validation, Trainings and Preventive Maintenance :

- HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.
- Trainings :We also take up preventive maintenance to reduce downtime of HPLC's Trainings.
- AMC's/CMC : We offer user training both in-House and at customer sites on PLC principles, operations, troubleshooting.
- Validations : We have protocols for carrying out periodic Validations as per GLP/ GMP/U SFDA norms.
- Instruments : We offer instruments/Renting Services Modules like pumps, detector etc. on Rent.



## ►► About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO :9001 Certified Company engaged in Designing, Manufaturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

## ►► Our Products & Technologies





## ►► Regulatory compliances



## ►► Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



# Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at [info@analyticalfoundation.org](mailto:info@analyticalfoundation.org)

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

## ►► Reach us @



# Analytical<sup>®</sup> Technologies Limited

**HPLC Solutions   MultipleLabs   Analytical Bio-Med   Analytical Distributors   Analytical Foundation (Trust)**

**Corporate & Regd. Office:**  
Analytical House, # E67 & E68,  
Ravi Park, Vasna Road, Baroda,  
Gujarat 390 015. INDIA

T: +91 265 2253620  
+91 265 2252839  
+91 265 2252370  
F: +91 265 2254395

E: [info@hplctechnologies.com](mailto:info@hplctechnologies.com)  
[info@multiplelabs.com](mailto:info@multiplelabs.com)  
[info@analyticalgroup.net](mailto:info@analyticalgroup.net)  
[info@analyticalbiomed.com](mailto:info@analyticalbiomed.com)

W: [www.ais-india.com](http://www.ais-india.com)  
[www.analyticalgroup.net](http://www.analyticalgroup.net)  
[www.hplctechnologies.com](http://www.hplctechnologies.com)  
[www.multiplelabs.com](http://www.multiplelabs.com)

Sales & Support Offices:  
across the country :  
Distributors & Channel  
partners World Wide

Note : Company reserves rights to add/delete/modify the contents /technical specifications of the catalogue without prior notice.