


Mobile X Ray Machine With Accessories

MXR 3029



EPC / PRODUCTS / APPLICATION / SOFTWARE / ACCESSORIES / CONSUMABLES / SERVICES

Analytical Technologies Limited

An ISO 9001 Certified Company

www.analyticalgroup.net

►► Specification

- High Frequency mobile X ray machine with output 100 mA or more. The mobile x ray equipment is required to perform x ray studies in emergency and trauma center and bedside in wards and ICU. The unit is compact, lightweight and easily transportable. It has following specifications.
- The unit is operative on mains voltage from single phase 180-260 v AC with automatic main compensation.
- Generator.
- Power : 4 kW
- kVp. Range : 40 – 100 kVp
m AS Range: 250 m As
- m A range : 10 mA to 100 mA
- Exposure Time: 10 ms to 5 sec.
- The digital display: kV and mAs parameters, System ON, System OFF, status and fault messages on the kV and mAs area.
- X RAY Tube: Stationary/Rotating Anode tube with focal spot 1.8 X 1.8 mm
- Tube stand: The tube stand is fully counterbalanced with rotation in all directions.
- Collimator: Collimator rotation is +90 to -90 degrees with auto shut off lamp facility.
- Cassette storage box: The equipment has cassette storage box for minimum of 4 cassette
- Ergonomics: The unit has small foot print. The height of the column stand is less than 150 cm for easy transportation in the lift etc. and areas with small height doors. The equipment is light weight, less than 130 kg.
- Brake system: The unit has effective brake system for parking.

►► HPLC Servicing, Validation, Trainings and Preventive Maintenance :

- HPLC Servicing : We have team of service engineers who can attend to any make of HPLC promptly @the most affordable cost.
- Trainings :We also take up preventive maintenance to reduce downtime of HPLC's Trainings.
- AMC's/CMC : We offer user training both in-House and at customer sites on PLC principles, operations, troubleshooting.
- Validations : We have protocols for carrying out periodic Validations as per GLP/ GMP/U SFDA norms.
- Instruments : We offer instruments/Renting Services Modules like pumps, detector etc. on Rent.



►► About Analytical Technologies

Analytical Technologies is synonymous for offering technologies for doing analysis and is the Fastest Growing Global Brand having presence in at least 96 countries across the globe. Analytical Technologies Limited is an ISO :9001 Certified Company engaged in Designing, Manufaturing, Marketing & providing Services for the Analytical, Chromatography, Spectroscopy, Bio Technology, Bio Medical, Clinical Diagnostics, Material Science & General Laboratory Instrumentation. Analytical Technologies, India has across the Country operations with at least 4 Regional Offices, 6 Branch Offices & Service Centers. Distributors & Channel partners worldwide.

►► Our Products & Technologies



►► Regulatory compliances



►► Corporate Social Responsibility

Analytical Foundation is a Nonprofit Organization (NGO) found for the purpose of:



Analytical Foundation

1. Research & Innovation Scientist's awards / QC Professional Award : Quality life is possible by innovation only and the innovation is possible by research only, hence ANALYTICAL FOUNDATION is committed to identify such personalities for their contributions across various field of Science and Technology and awarding them yearly. To participate for award, send us your details of research / testing / publication at info@analyticalfoundation.org

2. Improving quality of life by offering YOGA Training courses, Work shops / Seminars etc.

3. ANALYTICAL FOUNDATION aims to DETOXYFY human minds, souls and body by means of Yoga, Meditation, Ayurveda, Health Care, Awards, Media, Events, Camps etc.

►► Reach us @



 **Analytical**[®]
Technologies Limited

HPLC Solutions MultipleLabs Analytical Bio-Med Analytical Distributors Analytical Foundation (Trust)

Corporate & Regd. Office:
Analytical House, # E67 & E68,
Ravi Park, Vasna Road, Baroda,
Gujarat 390 015. INDIA

T: +91 265 2253620
+91 265 2252839
+91 265 2252370
F: +91 265 2254395

E: info@hplctechnologies.com
info@multiplelabs.com
info@analyticalgroup.net
info@analyticalbiomed.com

W: www.ais-india.com
www.analyticalgroup.net
www.hplctechnologies.com
www.multiplelabs.com

Sales & Support Offices:
across the country :
Distributors & Channel
partners World Wide

Note : Company reserves rights to add/delete/modify the contents /technical specifications of the catalogue without prior notice.